CrystalView™

ENGLISH


TABLE OF CONTENTS

Introduction: presentation of the CrystalView™ Irrigation Console	4
General description	5
Precautions for use	6
1.1 Operator profile	6
1.2 Target population	6
1.3 Intended use and contraindications	6
1.4 Side effects	6
1.5 Symbols Description	7
1.6 Warnings and general precautions	10
1.7 Technical characteristics	11
Product description	12
2.1 System indicators	12
2.2 Accessories and AC power cord	14
2.3 Consumable tubing	16
2.4 Audible indicators	18
2.5 Light indicators	18
2.6 Error messages	18
Receiving the equipment	19
3.1 Transport	19
3.2 Setup	19
3.3 System verification	19
Setup Procedures	20
4.1 Installing the console	20
4.2 Connecting the CrystalView™ Foot Pedal, 2-Way or 4-Way	20
4.3 Connecting the shaver interfaces	20
4.4 Connecting the CrystalView™ Interface Cables	22
4.5 Connecting the tubing	26
Tubing life cycle	31
5.1 After surgery	31
5.2 Connection of a new CrystalView™ Intermediary Tube Set for the next operation.	32
5.3 End of operating day	32
Console functions	33
6.1 Starting the console	33
6.2 Pressure level adjustment	33
6.3 Filling the pressure chamber	35
6.4 "Suction" mode	35
6.5 "Lavage" mode	36
6.6 Using the pedal	36
6.7 Using the earth contact	36

TABLE OF CONTENTS

Maintaining the device	37
7.1 Cleaning	37
7.2 Inspections (annual)	37
 Troubleshooting	 38
8.1 FAQ	38
8.2 Who to contact?	40
8.3 Distributor	40
8.4 Aftersales service and warranty	40
 Technical characteristics	 41
9.1 Physical dimensions	41
9.2 Storage conditions	41
9.3 Performance specifications	41
9.4 References	42
9.5 EMC compliance matrix	43

INTRODUCTION: PRESENTATION OF THE CRYSTALVIEW™ IRRIGATION CONSOLE

Dear Customer,

Thank you for purchasing your **CrystalView™ Irrigation Console**. This device is delivered with the technical documentation. Please always ensure that this manual is close hand; it describes your equipment in addition to its operation.

CrystalView™ is an irrigation console that can be used with its own tubing. Its system manages irrigation automatically as a function of suction, thereby providing a flow rate and a precisely controlled pressure. Its user-friendly and easy to use Human / Machine Interface (HMI) allows optimum control of the device for the practitioner. **CrystalView™ Irrigation Console** is intended for use during arthroscopy operations by orthopaedic surgeons and operating theatre nurses. The console was developed in collaboration with surgeons in order to meet their expectations and be capable of offering them maximum convenience.

The performance of the console allows:

- conveyance of the irrigation liquid from the bags of physiological saline to the joint via an arthroscopy sleeve, precisely controlling the pressure and flow rates,
- recovery of the contaminated fluids from the joint via a cannula or a surgical instrument to the waste collection bag,
- precise control of intra-articular pressure, regardless of the outgoing flow rate,
- selection of the appropriate suction level (proposed rates: Min, Low, Med, High),
- communication with the shaver consoles.

GENERAL DESCRIPTION

The **CrystalView™** system consists of a console, accessories and consumables.


Accessories:

- **CrystalView™** Foot Pedal, 2-Way and 4-Way
- **CrystalView™** Interface Clip
- AC Power cord
- **CrystalView™** Interface Cables

Consumable tubing:

- **CrystalView™** Day Use Tube Set
- **CrystalView™** Patient Tube Set
- **CrystalView™** Intermediary Tube Set

- A** **CrystalView™** Day Use Tube Set
- B** **CrystalView™** Patient Tube Set
- C** **CrystalView™** Intermediary Tube Set
- D** Arthroscope ConMed
- E** Shaver ConMed
- F** **CrystalView™** Foot Pedal
- G** Shaver Console ConMed
- H** Blue Cannula
- I** Shaver cable

Reading this manual is compulsory before any use of **CrystalView™ Irrigation Console**, its accessories and its consumable tubings.

PRECAUTIONS FOR USE

1

1.1 Operator profile

The **CrystalView™ Irrigation Console** is to be used in the operating theatre; the users include nurses and orthopedic surgeons.

1.2 Target population

Target population: adults and pediatric patients of more than 8 years of age, regardless of gender.

1.3 Intended use and contraindications

Indications

The **CrystalView™ Irrigation Console** is an arthroscopic pump system intended to provide fluid distension and irrigation of the knee, shoulder, hip, elbow, and ankle and wrist joint cavities and fluid suction during arthroscopic procedures.

Contraindications

The use of the **CrystalView™ Irrigation Console** is prohibited whenever arthroscopy is contraindicated.


1.4 Side effects

The known and possible side effects during arthroscopy procedure with a fluid management system are following :

Emphysema, hematoma, reflex sympathetic dystrophy, deep venous thrombosis (DVT) and fluid extravasations.

1.5 Symbols Description

Symbols for consoles and accessories

Symbol	Name	Description
	Fuse	External fuse rating: T5AH-250V Internal fuse rating: F2AL-250V
	Manufacturer	Identifies the manufacturer of the medical device.
	CE marking	Indicates that this equipment complies with the European medical device regulations. If applicable, there is also the identification of the notified body (four digits).
IP 20	Protection rating of the pump	2 = protected against solid foreign bodies of 12.5mm diameter and larger. 0 = no protection against water.
IP X7	Protection rating of the pedals	Protected against the effects of immersion (up to 1m).
	Medical device	Indicates that the product is a medical device.
	Quantity	Indicates the quantity of product present in the packaging.
	Temperature limit	Indicates the temperature limits to which the medical device can be safely exposed
	Humidity limit	Indicates the humidity limits to which the medical device can be safely exposed
	Distributor	Indicates the entity that distributes the medical device in the region concerned.
	Authorized representative in Switzerland	Specify the authorized representative in Switzerland.

PRECAUTIONS FOR USE

1

	Equipotential bonding symbol	Indicates presence of an earth contact on the rear of the CrystalView™ .
	Symbol for electrical and electronic equipment waste.	This device is subject to separate collection. Return the device to Hemodia at the end of its service life.
	Refer to the instruction manual	Reading the manual is compulsory before using the pump for the first time.
	Catalog reference	Indicates the manufacturer's catalog reference to allow formal identification of the medical device.
	Serial number	Indicates the serial number assigned by the manufacturer in order to formally identify a specific medical device.
	USB symbol	Indicates presence of a USB port on the rear of the CrystalView™ .
	Alternating current	Indicates that the CrystalView™ is to be connected to an alternating current power supply.
	Instruction manual; instructions for use	Informs the user to refer to the manual to obtain information about the CrystalView™ .
	Manufacturing batch	Indicates the number of the manufacturing batch.

PRECAUTIONS FOR USE

1

English

Tubing

	Use-by date	Indicates the date after which the medical device may no longer be used.
	Do not reuse	Indicates that a medical device may only be used once or on a single patient for one treatment only.
	Manufacturing date	Indicates the date on which the medical device was manufactured.
	Sterilized using ethylene oxide	Indicates that the medical device was sterilized with ethylene oxide.
	Manufacturing batch	Indicates the number of the manufacturing batch.

PRECAUTIONS FOR USE

1

1.6 Warnings and general precautions

The hospital staff are urged to read this manual before using and cleaning this product and its accessories. Failure to abide by these instructions may cause injuries and possible damages to the patient, or malfunction of the equipment.

The manufacturer of this device cannot be held liable for direct or consequential injuries or damage resulting from inappropriate use of the single-use products other than the Hemodia products. Any modification of the device, any repairs performed by an unauthorized service center or any use of single-use products other than the Hemodia products shall invalidate the manufacturer warranty and the civil liability coverage in case of material damages or bodily harm.

The equipment must only be used by personnel trained in arthroscopic procedures.

Do not connect the device to an unearthed or poorly earthed source of power. In order to avoid any risk of electric shock, this equipment must only be connected to a power supply equipped with a protective earth.

In order to avoid any risk of fire, replace spent fuses with fuses of the same characteristics

The system may be affected by electromagnetic interferences generated by other instruments. Check that the other instruments used in the operating theater comply with electromagnetic compatibility (EMC) standard IEC 60601-1-2.

If the system continues to suffer interferences, isolate the system and connect it to a different connector.

CrystalView™ Irrigation Console is a class 1 device in accordance with electrical safety standard NF EN 60601-1, the currently applicable version.

The pressure adjustments (chapter 6.2) suggested by ConMed are based on observations; the values in the table may be modified and should be adapted depending on the operation.

Risks of electrocution: do not remove the casing. Repairs may only be performed by qualified personnel.

Any serious incidents that have occurred in connection with these devices must be notified to the manufacturer and the competent authority of the Member State in which the user is based.

1.7 Technical characteristics

A - Sterilization

The **CrystalView™** Irrigation Console and the accessories (pedals, power cord, Shaver interface...) are not sterile. On the other hand, they must be decontaminated after each day of use.

On receiving the tubing, check that the packaging items have not been damaged. If a consignment is opened or damaged, or the expiry date is exceeded, contact the ConMed "Customer service" department by telephone or email.

This symbol applied to a packaging item means that the suction tubes are for single use. They cannot therefore be reused.

Retreatment is liable to change the characteristics of the materials, particularly by deforming and degrading the latter, which may have repercussions on the durability of the device and compromise its performances. These risks may potentially jeopardize patients' safety.

Never use a product if its sterile packaging is damaged or if the expiry date is exceeded.

B - Electrical characteristics

Input voltage: 100-240 VAC

Frequency: 50-60 Hz

Current draw: 500 VA

External fuses: T5AH-250V (2X)

Internal fuses: F2AL-250V (2X)

Overvoltage category: 2


Pollution level: 2

PRODUCT DESCRIPTION

2

2.1 System indicators

A - Front panel


- ❶ Transparent protective cover for the **CrystalView™** Irrigation Console.
- ❷ Irrigation pump head
- ❸ Pressure arm. Secures the tube around the console head.
- ❹ LCD screen. Displays the dynamic pressure, setpoint pressure and error messages.
- ❺ **ERROR.** The red indicator light illuminates when the device detects a problem. (refer to 2.6 Error messages).
- ❻ Button for selecting the pressure unit (mmHg and LH₂O).
- ❼ Transparent protective cover for the suction console.
- ❽ Suction console head.
- ❾ Retaining bracket for the suction tube.
- ❿ Indicator light for "lavage" mode.

- 11 The "Pinch Valve" trips automatically when the "shaver" is activated by selecting the appropriate suction tube.
- 12 **"Pinch Valve"** button. Opens the "Pinch Valve" automatically and pinches the cannula tube.
- 13 "On O/I" switch. Switches the device on or off. If the console is activated, "Console" appears on the display screen.
- 14 "Shaver" suction adjustment. Allows adjustment of the flow rate by the "shaver" (MIN, LOW, MED, HIGH).
- 15 **"Run/Stop"** button. Starts or stops the console. The green indicator light flashes when the console is stopped (set to stop).
- 16 Adjustment of suction by the cannula. Allows adjustment of the flow rate of the cannula (MIN, LOW, MED, HIGH).
- 17 Yellow connector allows to connect the **CrystalView™** Interface Clip.
- 18 Red connector allows connection of **CrystalView™ Foot Pedal, 2-Way or 4-Way**.
- 19 CPC™ connector. Allows connection of the **CrystalView™** Day Use Tube Set to the pressure sensor incorporated in the **CrystalView™ Irrigation Console**.
- 20 **"Fill Chamber"** button. Allows manual filling of the chamber and cancels the Low Pressure error.
- 21 Pressure adjustment button. 2 pressure units available:
In "LH₂O": press the arrows to increase or decrease the pressure level by increments/decrements of 5 LH₂O. Initial setting: 30 LH₂O Minimum level: 10 LH₂O / Maximum level: 120 LH₂O
In "mmHg": press the arrows to increase or decrease the pressure level by increments/decrements of 10 mmHg. Initial setting: 60 mmHg. Minimum level: 20 mmHg / Maximum level: 280 mmHg.
- 22 Retaining bracket for the irrigation tube.
- 23 Pressure chamber support.
- 24 Integrated "Shaver" interface. Detects functioning of the concurrent "shaver".
- 25 Blue connector (**CrystalView™** Interface Cables) allows interface through shaver operation through the **CrystalView™ Foot Pedal, 4-Way**.
- 26 Button to remove manually the tubing from "Pinch Valve".

PRODUCT DESCRIPTION

2

B - Back panel

- 27 Fuse holder and main socket
- 28 Equipotential terminal
- 29 Product reference number (REF)
- 30 Device serial number (SN)
- 31 Connector for software programming

2.2 Accessories and AC power cord

 Only connect the accessories dedicated for the **CrystalView™** Irrigation Console. Please, check the accessories available on page 42. Perform functional tests before use on the patient.

A - CV2044 or CV2042 CrystalView™ Foot Pedal

The 4- or 2-way pedals are connected to the **CrystalView™** Irrigation Console. They allow activation of suction by the "cannula" and "shaver" tubes at the preconfigured levels and activation or deactivation of "lavage" mode.

- 32 The black "**LEFT**" pedal makes it possible to control a "shaver" from a shaver console when the **CrystalView™** Interface Cable is connected and activate suction by the "shaver" at a preconfigured level.
- 33 The black "**RIGHT**" pedal makes it possible to control a "shaver" from a shaver console when the **CrystalView™** Interface Cable is connected and activate suction by the "shaver" at a preconfigured level.
- 34 Blue "**Cannula**" pedal. Activates suction by the cannula at the preconfigured level or cancels "lavage" mode.
- 35 Red "**Lavage**" pedal. Activates "lavage" mode.

B - Interfaces (refer to § 4.3 and 4.4)

There are three different types of interfaces:

- The **CrystalView™** Interface Clip: detects operation of a "shaver" from a console and activates "Shaver" suction at the preconfigured level
- The Integrated "Shaver" interface (incorporated in the front panel): detects operation of a "shaver" from a console and activates "Shaver" suction at the preconfigured level.
- The **CrystalView™** Interface Cables: controls the shaver hand pieces and activates "Shaver" suction at the preconfigured level from 4 way footpedal.

C - AC power cord

The AC power cord is used to connect the **CrystalView™** Irrigation Console to the power supply.

Please use the appropriated AC power cord to any local AC mains outlet.


The AC power cord can be used as a means of emergency stop; disconnecting the AC power cord makes it possible to isolate the **CrystalView™** Irrigation Console from the power supply.

Install the **CrystalView™** Irrigation Console such that it is easy to disconnect the AC power cord from the power supply.

2.3 Consumable tubing


CrystalView™ Day Use Tube Set (Ref: CV2025)

This tubing connects the irrigation bags to the **CrystalView™** Intermediary Tube Set. The pressure is read through a filter connected to the **CrystalView™** Irrigation Console. The **CrystalView™** Day Use Tube Set irrigation is a tubing for daily use. It can be used for several consecutive operations for different patients. The duration of use must not however exceed 14 hours. No contamination is possible among patients if, and only if, the tubing is used in conjunction with an intermediary tubing equipped with a check valve. The procedure for installing and changing the tubing is explained in chapters 4.5 and 5.

CrystalView™ Intermediary Tube (Ref: CV2020)


The **CrystalView™** Intermediary Tube is used to convey the sterile saline solution from the **CrystalView™** Day Use Tube to the arthroscope.

 The tubing is intended for single use and must not be reused or re-sterilized.

CrystalView™ Patient Tube Set (Ref : CV2030)

This tubing allows suction from the joint, either via a cannula or via another suction device "shaver". Both modes are controlled by the "Pinch Valve".

 The tubing is intended for single use and must not be reused or re-sterilized.

PRODUCT DESCRIPTION

2

2.4 Audible indicators

A double signal sounds when the machine is switched on (when the "On O/I" button is pressed) and when "lavage" mode is stopped.

A single audible signal sounds each time a button is pressed.

2.5 Light indicators

- 5 This red indicator light illuminates when a safety parameter is not observed.
- 6 This illuminated indicator light shows the pressure unit selected ("mmHg" or "LH₂O").
- 10 This indicator light illuminates when "lavage" mode is active.
- 12 This indicator light illuminates when the "Pinch Valve" is active.
- 14 16 These illuminated indicator lights show the flow rates selected (MIN, LOW, MED, HIGH).
- 15 This indicator light illuminates when the console is set to "Run" mode; the indicator light flashes when the console is set to "Stop" mode.
- 20 This indicator light illuminates when the pinch valve is active.

2.6 Error messages

Displayed on the screen	Description
LOW PRESSURE	Dynamic pressure is less than 10 mmHg (or 5 LH ₂ O).
Leak of Pressure	A leak has been detected in the irrigation tube.
Sensor Pressure Error	Pressure measurement error.

3.1 Transport

If the device is dropped or in case of damage to the equipment and its accessories, do not use the console or any of its accessories and contact the Hemodia Aftersales Service by telephone or email (maintenance@hemodia.com or 00 44 20 34 45 51 79).

3.2 Setup

- Connect the AC power cord supplied by ConMed to the mains socket located at the back of the console. 27
- Ensure that the cable is properly inserted into its socket.
- Connect the shaver of the console depending on the selected mode according to the connections indicated on chapter "4.3 et 4.4".
- Connect the pedal supplied by ConMed to the connector of the console. 18

3.3 System verification

- Connect the AC power cord to a mains supply.
- Switch on the console by setting the "**On I/O**" switch to the "I" position. 13
- Start the console by pressing the "**Run/Stop**" button 15, the message "**PUMP**" appears on the LCD screen. 4
- Press all the buttons (each activated button should emit an audible signal), with the exception of the "On" button, 13 which emits 2 audible signals on start-up.

SETUP PROCEDURES

4

4.1 Installing the console

 Install the console on a level surface to avoid any risk of the machine falling.
To avoid any risk of interferences with the integrated shaver interface, do not place the **CrystalView™** Irrigation Console directly on the RF console.
Connect the power supply cable supplied by ConMed to an electrical socket.
Raise both covers (irrigation and suction).
Press the "On O/I" switch and set it to "I". The LCD screen illuminates and displays "Console".

4.2 Connecting the CrystalView™ Foot Pedal, 2-Way or 4-Way

Connect the **CrystalView™** Foot Pedal, 2-Way or 4-Way supplied by ConMed to the red connector of the console. **18**

4.3 Connecting the shaver interfaces

A - CrystalView™ Interface Clip


Connect the **CrystalView™** Interface Clip to the "Shaver" connector of the console. **17**
Connect the clip of the **CrystalView™** Interface Clip around the handpiece cable.
Ensure that the clip is as close as possible to the "Shaver" console.
Connect the pedal to the pedal connector of the "Shaver" console.

NOTE : In the case where the Shaver clip interface malfunctions, you can use the integrated Shaver.

The Shaver clip interface is compatible with the following shaver console/handpiece systems:


- * **CONMED/LINVATEC** - Console D3000 / D4000 - Shaver ERGO / ADVANTAGE TURBO
- * **SMITH & NEPHEW/DYONICS** - Console EP1, POWER, POWER II - Shaver POWER-MAX ELITE / POWER / POWERMAX
- * **STRYKER** - Console CORE - Shaver FORMULA 180 / FORMULA / 12K SHAVER
- * **ARTHREX** - Console APS II - Shaver AR-8325H / AR-8325F / AR-8320F / AR-8330F

! Do not use the **CrystalView™** Interface Clip and the **CrystalView™** Interface Cable at the same time.


B - Integrated shaver interface

Place the handpiece cable in the fixing bracket **24**.
Connect the Shaver pedal to the pedal connector of the Shaver console.

! Do not use the integrated Shaver interface and **CrystalView™** Interface Clip at the same time.

4.4 Connecting the CrystalView™ Interface Cables

Console and Handpiece activation

INSTRUCTIONS FOR USE

Connect fully the **CrystalView™** Interface Cable into the Blue connector of the **CrystalView™** Irrigation Console. **25**

Then, from **CrystalView™** Irrigation Console, connect the **CrystalView™** Interface Cable into the foot connector of the adapted shaver console. **A**

PRECAUTIONS

Before plugging and unplugging the **CrystalView™** Interface Cable, turn off both electromedical devices. Prior to use, inspect the product and do not use if damaged. Control the **CrystalView™** Interface Cable connectors and ensure that cleaned cable connectors are completely dry prior to connecting the electromedical devices.

Do not use the **CrystalView™** Interface Cable and the integrated Shaver interface at the same time.

Do not use the **CrystalView™** Interface Cable and the **CrystalView™** Interface Clip at the same time.



A - CV2050 (CrystalView™ Interface Cable for ConMed shaver)

The **CrystalView™** Interface Cable for ConMed shaver (REF: CV2050) is used during arthroscopy procedure to activate a **ConMed / Linvatec** handpiece of a **ConMed / Linvatec** shaver console using a **CrystalView™** Foot Pedal, 4-Way. The **CV2050** is compatible with following shavers consoles : **D3000 ADVANTAGE/ D4000 from ConMed/ Linvatec**.

 It is recommended before using the **CrystalView™** Interface Cable for ConMed shaver, to read the instructions ConMed/Linvatec shaver console user manual for the warnings and precautions associated with those devices.

When using the **ConMed / Linvatec** shaver console, the activation of the console and handpiece can be made from the **CrystalView™** Foot Pedal, 4-Way:

Press the **LEFT** from **CrystalView™** Foot Pedal, 4-way to choose rotation mode (FORWARD / REVERSE / OSCILLATION) on ConMed / Linvatec shaver console.

 mode FORWARD	 mode REVERSE	 mode OSCILLATION
---	---	---

Press and keep pressed **RIGHT** from **CrystalView™** Foot Pedal, 4-way to activate the Handpiece

 For each operating mode (**FORWARD/REVERSE/OSCILLATION**), verify that the system is functioning properly before using it outside the joint.

SETUP PROCEDURES

B - CV20SN (CrystalView™ Interface Cable for Smith and Nephew shaver)

The **CrystalView™** Interface Cable for Smith and Nephew shaver (REF: CV20SN) is used during arthroscopy procedure to activate a SMITH & NEPHEW / DYONICS handpiece of a SMITH & NEPHEW / DYONICS shaver console by using the **CrystalView™** Foot Pedal, 4-Way.

The reference CV20SN is compatible with the following shavers consoles: EP1 / POWER / POWER II Control System from SMITH & NEPHEW / DYONICS.

It is recommended before using the the **CrystalView™** Interface Cable for Smith and Nephew shaver to read the instructions SMITH & NEPHEW / DYONICS shaver console user manual for the warnings and precautions associated with those devices.

When using the SMITH & NEPHEW / DYONICS shaver console, the activation of the pump and handpiece can be made from the **CrystalView™** Foot Pedal, 4-Way::

Press and keep pressed RIGHT from CrystalView™ Foot Pedal, 4-Way to activate the handpiece in Forward mode.	 Forward mode
Press and keep pressed LEFT from CrystalView™ Foot Pedal, 4-Way to activate the handpiece in Reverse mode.	 Reverse mode
Press and keep pressed LEFT and RIGHT from CrystalView™ Foot Pedal, 4-Way to activate the handpiece in the Oscillation mode.	 Oscillation mode

For each operating mode (FORWARD/REVERSE/OSCILLATE), verify that the system is functioning properly before using it outside the joint.

C - CV20AR (CrystalView™ Interface Cable for Arthrex shaver)

The **CrystalView™** Interface Cable for Arthrex shaver (REF: CV20AR) is used during arthroscopy procedure to activate an ARTHREX handpiece of an ARTHREX shaver console by using the **CrystalView™** Foot Pedal, 4-Way. The reference CV20AR is compatible with following shavers consoles: POWER SYSTEM II / SYNERGY from ARTHREX.

 It is recommended before using the **CrystalView™** Interface Cable for Arthrex shaver to read the instructions ARTHREX shaver console user manual for the warnings and precautions associated with those devices.

When using the ARTHREX shaver console, the activation of the pump and handpiece can be made from the **CrystalView™** Foot Pedal, 4-Way.

Press and keep pressed RIGHT from CrystalView™ Foot Pedal, 4-Way. to activate the handpiece in Forward mode.	 Forward mode
Press and keep pressed LEFT from CrystalView™ Foot Pedal, 4-Way. to activate the handpiece in Reverse mode.	 Reverse mode
Press and keep pressed LEFT and RIGHT from CrystalView™ Foot Pedal, 4-Way to activate the handpiece in the Oscillation mode.	 Oscillation mode

 For each operating mode (FORWARD/REVERSE/OSCILLATE), verify that the system is functioning properly before using it outside the joint.


SETUP PROCEDURES

4

4.5 Connecting the tubing

A - CrystalView™ Day Use Tube Set (REF: CV2025)

Indications: This tubing for daily use delivers the saline solution (physiological saline) from the bag to the **CrystalView™** Intermediary Tube Set. This device is used in conjunction with the **CrystalView™** Intermediary Tube Set and the **CrystalView™** Irrigation Console.



Method of use: Before any handling operations, close the three white clamps and connect the spikes to the bags of saline solution.




Slip the green coupling of the **CrystalView™** Day Use Tube Set into the mounting bracket provided for the irrigation tubing.

22

Position and center the tube around the irrigation console. Pull the tubing tight. Place the pressure chamber (vertically) in its support located on the machine. The pressure chamber must remain vertical during the operation and it is recommended to position that the console on the same level as the patient.

Close the pressure arm on the console body in order to keep the tubing around the console head

Press the CPC™ connector button of the console. Push the hydrophobic filter to connect properly the tubing of the Day Set into the CPC™ connector of the console. 19


Open one of white clamp under spike. Close the left-hand protective cover 1 and press the «Run/Stop» button. 15


Press the "Fill Chamber" button several times 20 in order to fill the pressure chamber up to the mark engraved on the pressure chamber (1/3 of its capacity).

B - CrystalView™ Intermediary Tube Set (REF: CV2020)

Indications: This tubing delivers the physiological saline from the day irrigation tube to the arthroscope.

Method of use:

Use an aseptic method. After having checked that the "Luer-Lock" coupling is correctly screwed in place, the scrub nurse hands the blue end to the circulating nurse.

SETUP PROCEDURES

4

The circulating nurse removes the blue caps from the **CrystalView™** Intermediary Tube Set and the **CrystalView™** Day Use Set and immediately connects the couplings of these two tubings. Firmly tighten both "Luer-Lock" couplings.


The console must be set to "Run" mode. Open the clamp under the pressure chamber in order to fill the irrigation tubing. Close off the clamp under the pressure chamber when the irrigation tubing is purged.

C - CrystalView™ Patient Tube Set (REF: CV2030)

The **CrystalView™** Patient Tube Set (CV2030) is composed with **CrystalView™** Intermediary Tube Set (CV2020) and a suction tubing.

Indications: this tubing is used to aspirate the fluids in the joint via the cannula or a shaver.

Method of use:

Use an aseptic method. The scrub nurse hands the tube end marked "waste" to the circulating nurse who inserts it into the waste collection system.

SETUP PROCEDURES





4

English

The circulating nurse slips the orange connector of the suction tubing into the orange mounting bracket of the console.


Position and center the tube around the console head.

Pull the tubing tight, insert it into the recess provided and close the pressure arm on the console body in order to keep the tube around the console head.

Insert the right-hand "outflow cannula" tube into the right side of the "Pinch Valve" 11 and push it to the back.

Press the "Pinch Valve" button and insert the "shaver" tube into the left side of the "Pinch Valve" 11 and push it to the back.

SETUP PROCEDURES

4

Pass the tubing exiting the clamp through the hook provided for this purpose.
Close the right-hand protective cover. **7**



Handling operations to be performed in the sterile area


Insert the cannula into the OUTFLOW CANNULA tube.

Plug the "shaver" tube on to the suction connector of the handpiece.

Connect the "Luer-Lock" coupling of the irrigation to the arthroscope before opening the stopcock of the arthroscope.

Open the white clamp under the pressure chamber on the irrigation side.

5.1 After surgery

Leave the **CrystalView™** Day Use Tube Set in place when another operation is scheduled on the same day.

Close the clamp under the pressure chamber.

Press the blue "Cannula" 34 pedal in order to drain the joint.

Disconnect the arthroscope.

Disconnect the white "Luer-Lock" coupling from the **CrystalView™** Intermediary Tube Set.

The check valve is to remain connected to the **CrystalView™** Day Use Tube Set until the next operation in order to ensure sterility between procedures.

Suction tubing withdrawal:


Close all the white clamps.

Press the blue "**Cannula**" pedal in order to drain the joint and the suction tubing.

Press the "**RIGHT**" or "**LEFT**" pedals in order to drain the "shaver" tubing.

Remove the "shaver" handpiece and the cannula. Remove the right-hand "outflow cannula" tubing located on the "Pinch Valve".

Press the "**Pinch Valve**" button or press the "**Manual Pinch Valve**" button to remove the "shaver" tubing.

The tubing is designed to be single use. Retreatment is liable to change the characteristics of the materials, particularly by deforming and degradation, which may have repercussions on the durability of the device and compromise its performances. These risks may potentially jeopardise patients' safety.



Refer to the waste disposal procedure applicable in your institution in order to dispose of the tubing.

TUBING LIFE CYCLE

5

5.2 Connection of a new CrystalView™ Intermediary Tube Set for the next operation.


After having checked that the "Luer-Lock" is correctly locked, the scrub nurse hands the blue end to the circulating nurse.

The circulating nurse unscrews the blue "Luer-Lock" coupling on the **CrystalView™** Day Use Tube Set irrigation tubing removes the check valve still in place and immediately connects the "Luer-Lock" coupling on the **CrystalView™** Day Use Tube Set to the blue "Luer-Lock" of the new intermediary irrigation tubing.

Never connect two check valves.


5.3 End of operating day

Close all the white clamps.

Press the blue "**Cannula**" pedal in order to drain the joint and the suction tubing.

Press the "**RIGHT**" or "**LEFT**" pedals in order to drain the "shaver" tubing of the suction tubing.

Press the CPC™ connector to remove the tube with the hydrophobic filter.

Remove all the tubing. Remove the "shaver" handpiece and the cannula.

Refer to the waste disposal procedure applicable in your institution in order to dispose of the tubing.

Switch off the console by setting the "**On I/O**" switch to the "O" position.

Clean the **CrystalView™** devices (Console and accessories) after each day of use.

6.1 Starting the console

Set the "On O/I" switch to "I" 13 and press the "Run/Stop" button. 15

6.2 Pressure level adjustment

A - Adjustment

The pressure is displayed on the LCD screen of the console 4. Press once on the pressure adjustment buttons (+ or -) 21 allows visualization of the preset pressure. Press a second time to adjust preset pressure.

The default pressure setting is 30 LH₂O.

It is recommended to start the lowest possible pressure setting, so as to obtain the desired expansion in the joint. The intra-articular pressure can then be increased.

Reminder:

In mmHg: press the arrows to increase or decrease the pressure level by increments/decrements of 10. Minimum level: 20 / Maximum level: 280.

In LH₂O: press the arrows to increase or decrease the pressure level by increments/decrements of 5. Minimum level: 10 / Maximum level: 120.

B - Adjustment of Cannula and Shaver flow rates

Four cannula suction levels (MIN, LOW, MED, HIGH) can be set on the **CrystalView™ Irrigation Console**. Press on the "**CANNULA SUCTION**" 16 button to adjust the cannula suction flow rate.

Four Shaver suction levels (MIN, LOW, MED, HIGH) can be set on the **CrystalView™ Irrigation Console**. Press on the "**SHAVER SUCTION**" 14 button to adjust the shaver suction flow rate.

C - Recommended pressure and flow rate level

Joint	Pressure level without tourniquet		Cannula suction	Shaver suction
	mmHg	LH ₂ O		
Shoulder	40	20	Low or Med	Low or Med
Acromioplasty	60	30	Low or Med	Med or High
Knee	30	10	Low or Med Low or Med	Low or Med Med or High
Wrist	20	10	Low or Med	Low or Med
Elbow, ankle	20	10	Low or Med	Low or Med
Hip	50	25	Low or Med	Low or Med

CONSOLE FUNCTIONS

D - Pressure unit conversion

mmHg	LH ₂ O	LH ₂ O	mmHg
20	8,9	10	22,5
30	13,3	15	33,8
40	17,8	20	45,0
50	22,2	25	56,3
60	26,7	30	67,5
70	31,1	35	78,8
80	35,6	40	90,0
90	40,0	45	101,3
100	44,4	50	112,5
110	48,9	55	123,8
120	53,3	60	135,0
130	57,8	65	146,3
140	62,2	70	157,5
150	66,7	75	168,8
160	71,1	80	180,0
170	75,6	85	191,3
180	80,0	90	202,5
190	84,4	95	213,8
200	88,9	100	225,0
210	93,3	105	236,3
220	97,8	110	247,5
230	102,2	115	258,8
240	106,7	120	270
250	111,1		
260	115,6		
270	120,0		
280	124,4		

E - Errors

E1 - High pressure error

The "**ERROR**" indicator light illuminates if the pressure increases beyond the safety level (140 LH₂O or 300 mmHg). The irrigation console stops until the setpoint pressure is restored.

E2 - Low pressure error

The "**ERROR**" indicator light illuminates if the pressure falls below the safety level (5 LH₂O or 10 mmHg), "**LOW PRESSURE**" is displayed on the LCD screen and the console stops.

Check that the CPC™ coupling is properly connected (if not reconnect the CPC™ coupling and press the "**Run/Stop**" button to restart the console).

Check if the hydrophobic filter is wet (in this case replace the **CrystalView™** Day Use Tube Set and press the "**Run/Stop**" button to restart the console).

E3 - Other types of error

When an error is detected, the error indicator light illuminates and the error code flashes on the LCD screen. Refer to table 2.6 on page 18.

6.3 Filling the pressure chamber

The "**Fill Chamber**" button 20 is used to fill the pressure chamber. The pressure chamber must be filled to one third of its capacity.

6.4 "Suction" mode

The "**Cannula Suction**" 16 and "**Shaver Suction**" buttons 14 allow adjustment of the suction flow rate of the cannula or "shaver". There are 4 settings: Min, Low, Med, and High.

An indicator light is illuminated opposite the current setting.

The "Pinch Valve" allows selection of the suction mode: via the cannula or "shaver".

A - Via the cannula

Suction mode via the cannula is activated by pressing the blue "**Cannula**" pedal 34; the indicator light corresponding to the selected preconfigured level flashes.

B - Via the shaver

Suction mode via the "shaver" is activated by pressing the "**LEFT**" pedal 32 or the "**RIGHT**" pedal 33 or both pedals at the same time. The indicator light corresponding to the selected preconfigured level flashes.

6.5 "Lavage" mode

"Lavage" mode is used to limit bleeding.

Press the red "**Lavage**" pedal **35** for 1 second to activate "lavage" mode. The green "Lavage" indicator light **10** remains illuminated throughout the lavage cycle.

The pressure increases to 50% from the preset pressure for 2 minutes.

Press the blue "Cannula" pedal to cancel a lavage cycle. **34**

An audible beep sound indicates activation of "lavage" mode; two audible beep sounds indicate that "lavage" mode is terminated.

6.6 Using the pedal

- Using the pedal: With the ConMed Shaver Ergo™ Handpieces and Advantage™ Turbo handpieces, the "**LEFT**" black pedal selects the mode (Forward / Reverse / Oscillate) and the "**RIGHT**" black pedal activates suction via the «shaver» and activates the handpiece in selected mode.
- Using the pedal: With the SMITH & NEPHEW / DYONICS, ARTHREX Shaver, the "**LEFT**" black pedal activates suction via the «shaver» and activates the handpiece in reverse mode. The "**RIGHT**" black pedal activates suction via the «shaver» and activates the handpiece in forward mode. Pressing the black "**RIGHT**" and "**LEFT**" pedals simultaneously activates suction via the «shaver» and activates the handpiece in oscillation mode.

6.7 Using the earth contact

The earth contact allows connection of the **CrystalView™ Irrigation Console** to another device via a connecting cable in order to equalize the potentials and avoid any risk to the patient in case of a fault in the electrical installations. The connecting cable is not supplied.

7.1 Cleaning

At the end of the operating day, switch off the **CrystalView™** Irrigation Console, disconnect the device and remove the used tubing. Clean all the surfaces of the console in addition to the protective covers using a cloth moistened with a neutral pH detergent. Wipe the device afterwards with a cloth soaked in distilled water. Clean all the surfaces of accessories (AC Power cord, Foot Pedals, interfaces) with a neutral pH detergent.

- Do not use flammable agents for cleaning or disinfection.
- Do not immerse the **CrystalView™** Irrigation Console, the interfaces, AC Power Cord and the electrical connectors of the pedals during cleaning.
- Observe the cleaning procedures.
- Do not autoclave the **CrystalView™** Irrigation Console and its accessories.
- Do not fold the pedal cable and interface when stowing away.

In order to avoid electric shock, disconnect the power supply before cleaning the console.

- The electrical connections must remain dry.
- Do not fold the mains cable during storage.

7.2 Inspections (annual)

Hemodia recommends inspecting your equipment once a year to check serviceability of the console in addition to its technical compliance. Only Hemodia may perform this annual inspection.

TROUBLESHOOTING

8.1 FAQ

Problem	Cause	Solution
No pressure available.	The pressure connector is not connected.	Connect the CPC™ pressure connector to the console.
The console stops. The red " ERROR " indicator light is illuminated.	Low pressure safety has triggered.	Check that the CPC™ coupling is properly connected. Press the " Run/Stop " button to restart the console.
	The hydrophobic filter is wet.	Replace the "Day Set" tubing.
	There is no more irrigation liquid.	Replace the bag of saline solution and press the " Fill Chamber " key.
No suction via the cannula.	The cannula is blocked.	Reposition the cannula or clean it.
	The red clamp is closed.	Open the red clamp.
No suction via the "shaver".	The "Pinch Valve" is faulty.	Check the positioning of the suction tubes. Check the interface connecting cable.
	The "shaver" blade is clogged.	Disassemble the blade and clean or change it.
The console heads do not rotate.	The " On/Off " switch is set to "O" or the console is set to "stop". The cover are opened.	Check that the " On/Off " switch is set to "I" and that the console is in "run" mode. Close both protective covers.
There is no irrigation.	The stop cock of the arthroscope or a clamp is closed, preventing irrigation.	Check the irrigation tube from the bags to the joint. Open the closed clamps or cocks.

Problem	Cause	Solution
The irrigation console rotates very rapidly and emits excessive noise.	No irrigation liquid circulation.	Replace the bags of saline solution and/or check the clamps under the bags are opened. Press the "Fill Chamber" button if the pressure chamber is less than 1/3 full.
The irrigation console surges, resulting in abnormal pressure fluctuations.	Water has maybe backed up into the tube as far as the filter.	If water has entered the pressure detection duct, change the CrystalView™ Day Use Tube Set.
Overfilling of the pressure chamber.	Air leak in the irrigation tube.	Check that the CPC™ connector is well connected or undamaged. Change the CrystalView™ Day Use Tube Set if necessary and restart the console.
The pressure in the system is inadequate.	An irrigation problem is present.	Check the arthroscope stop cock, the clamps under the bags of saline solution and under the pressure chamber are opened.
The suction in the system is inadequate.	The cannula has been poorly positioned in the subcutaneous tissues or is obstructed.	Reposition the cannula or clean it. Check that the arthroscope sheath provides an adequate flow rate.
The CrystalView™ Irrigation Console remains stuck in Shaver suction.	The electronic card or Pinch Valve is defective.	Turn off the CrystalView™ Irrigation Console so you can remove the CrystalView™ Patient Tube Set (CV2030). Contact Aftersales Service.

8.2 Who to contact?

Hemodia is the legal manufacturer. Please contact the Hemodia Aftersales Service by email or telephone if a problem arises.

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8.4 Aftersales service and warranty

The warranty period for the **CrystalView™** Irrigation Console is 12 months. Within this period, errors resulting from faulty material and/or inadequate workmanship will be remedied by the manufacturer free of charge.

Transport costs and shipping risks, however, will not be carried by the manufacturer.

The manufacturer is not liable for direct or consequential damages, and the warranty becomes null and void if:

- the device and/or the accessories are improperly used, prepared, or maintained;
- the instructions and rules in the instructions for use are not adhered to;
- unauthorized persons perform repairs, adjustments, or alterations on the device or accessories;
- unauthorized persons open the device;
- the prescribed inspection and maintenance schedule is not adhered to.

Only certified service providers are allowed to perform repairs or alterations to the device or accessories.

9.1 Physical dimensions

Dimensions

Height: 202 mm

Width: 434 mm

Depth: 377 mm

Weight: 11 kg

9.2 Storage conditions

Temperature range: between -10°C and +50°C.

Relative humidity range: included between 30% and 90%.

9.3 Performance specifications

- Pressure level:

In mmHg: minimum level: 20 / maximum level: 280

In LH₂O: minimum level: 10 / maximum level: 120

- Inflow rate:

Min: 100 ml/min. (+/- 50 ml/min.)

Max: 1200 ml/min. (+/- 100 ml/min.)

- Cannula flow rate:

Resting: 100 ml/min. +/- 50 ml/min.

Min: 200 ml/min. +/- 50 ml/min

Low: 300 ml/min. +/- 75 ml/min.

Med: 500 ml/min. +/- 75 ml/min.

High: 600 ml/min. +/- 75 ml/min.

- Shaver flow rate:

Min: 200 ml/min. +/- 50 ml/min.

Low: 400 ml/min. +/- 75 ml/min.

Med: 600 ml/min. +/- 75 ml/min.

High: 800 ml/min. +/- 100 ml/min.

- Software Version

X10.12

9.4 References

Product ref.	Product description
CV2000	CrystalView™ Irrigation console
CV2044	CrystalView™ Foot Pedal, 4-Way
CV2042	CrystalView™ Foot Pedal, 2-Way
CV2051	CrystalView™ Interface Clip
CV2025	CrystalView™ Day Use Tube Set
CV2020	CrystalView™ Intermediary Tube Set
CV2030	CrystalView™ Patient Tube Set
103773	AC Power Cord for EU
103822	AC Power Cord for UK
103899	AC Power Cord for CH
CV2050	CrystalView™ Interface Cable for ConMed shaver
CV20SN	CrystalView™ Interface Cable for Smith and Nephew shaver
CV20AR	CrystalView™ Interface Cable for Arthrex shaver
15801580	Arthroscopic cannula 4.5mm (Blue Cannula)

9.5 EMC compliance matrix

Electromagnetic immunity - Enclosure port, AC input power port, DC input power port, patient coupling port and input/output signal ports Professional healthcare facility environment			
The CrystalView™ pump is designed for use in the electromagnetic environment specified below. The client or user of the CrystalView™ pump should ensure that it is used in such an environment.			
Basic EMC phenomenon/standard or test method	Test levels	Compliance level	Type of port
Electrostatic discharges (ESD)/IEC 61000-4-2	$\pm 8 \text{ kV} \rightarrow \text{contact}$ $\pm 2 \text{ kV}, \pm 4 \text{ kV}, \pm 6 \text{ kV},$ $\pm 8 \text{ kV}, \pm 15 \text{ kV} \rightarrow \text{air}$	Indirect $\pm 8 \text{ kV}$ (vertical) Indirect $\pm 8 \text{ kV}$ (horizontal) Direct $\pm 8 \text{ kV}$ Air $\pm 2 \text{ kV}$ Air $\pm 4 \text{ kV}$ Air $\pm 8 \text{ kV}$ Air $\pm 15 \text{ kV}$	- Enclosure port - Patient coupling port - Input/output signal ports
Radiated RF electromagnetic fields/IEC 61000-4-3	3 V/m 80 MHz - 2.7 GHz 80% AM at 1 kHz	80 MHz - 1 GHz $\rightarrow 3 \text{ V/m}$ (80% AM, 1kHz) on 4 sides in vertical and horizontal polarity 1 GHz - 2.7 GHz $\rightarrow 3 \text{ V/m}$ (80% AM, 1kHz) on 4 sides in vertical and horizontal polarity	- Enclosure port
Proximity fields emitted by wireless RF communication devices/ IEC 61000-4-3	9 V/m 704 MHz - 787 MHz and 5.1 GHz - 5.8 GHz PM ext. 217 Hz 27 V/m 380 MHz - 390 MHz PM ext. 18 Hz 28 V/m 430 MHz - 470 MHz FM deviation +5 kHz sinus. 1 kHz 28 V/m 860 MHz - 960 MHz PM ext. 18 Hz 28 V/m 1.7 GHz - 1.99 GHz and 2.4 GHz - 2.57 GHz PM ext. 217 Hz	9 V/m 710 MHz, 745 MHz, 780 MHz and 5,24 GHz, 5,5 GHz, 5,785 GHz 27 V/m 385 MHz 28 V/m 450 MHz 28 V/m 810 MHz, 870 MHz, 930 MHz 28 V/m 1,72 GHz, 1,845 GHz, 1,97 GHz and 2,45 GHz	- Enclosure port
Surges between phase and neutral/IEC 61000-4-5	$\pm 0.5 \text{ kV}, \pm 1 \text{ kV}$	AC power supply with dephasing of $0^\circ, 90^\circ, 180^\circ, 270^\circ$ $\pm 0.5 \text{ kV}$ (level 2) $\pm 1 \text{ kV}$ (level 3)	- AC input power port - DC input power port

TECHNICAL CHARACTERISTICS

Electromagnetic immunity - Enclosure port, AC input power port, DC input power port, patient coupling port and input/output signal ports Professional healthcare facility environment			
The CrystalView™ pump is designed for use in the electromagnetic environment specified below. The client or user of the CrystalView™ pump should ensure that it is used in such an environment.			
Basic EMC phenomenon/standard or test method	Test levels	Compliance level	Type of port
Surges between phase and earth and then between neutral and earth/IEC 61000-4-5	$\pm 0.5 \text{ kV}$, $\pm 1 \text{ kV}$, $\pm 2 \text{ kV}$	AC power supply with dephasing of 0° , 90° , 180° , 270° $\pm 0.5 \text{ kV}$ (level 1) $\pm 1 \text{ kV}$ (level 2) $\pm 2 \text{ kV}$ (level 3)	- AC input power port - DC input power port - Input/output signal ports
Conducted disturbances, induced by RF fields/IEC 61000-4-6	3 V 0.15 MHz - 80 MHz 6 V in ISM bands and bands between 0.15 MHz and 80 MHz 80% AM at 1 kHz	AC power supply 230Vac/50Hz Pedal cable (Pedal side) Pedal cable (Pump side) 3V on the band between 0.15 MHz and 80 MHz 6V on the frequency bands: 6.765 MHz - 6.795 MHz; 13.553 MHz - 13.567 MHz; 26.957 MHz - 27.283 MHz; 40.66 MHz - 40.70 MHz.	- AC input power port - DC input power port - Input/output signal ports - Patient coupling port
Magnetic fields at the designated industrial frequency/IEC 61000-4-8	30 A/m 50 Hz or 60 Hz	Not applicable No magnetically sensitive components or circuits	- Enclosure port
Voltage dips/IEC 61000-4-11	0% UT; 0.5 cycle At 0° , 45° , 90° , 135° , 180° , 225° , 270° and 315° 0% UT; 1 single-phase cycle: at 0° 70% UT; 25/30 single-phase cycle: at 0°	AC power supply at: - 0% UT (0° , 45° , 90° , 135° , 180° , 225° , 270° and 315°) for 10 ms - 0% UT (0°) for 20 ms - 70% UT (0°) for 500 ms	- AC input power port
Voltage interruptions/IEC 61000-4-11	0 % UT; 250/300 cycle	AC power supply at 0% UT (0°) for 5 sec.	- AC input power port
Transient electrical disturbances by conduction along supply lines only/ISO 7637-2	Not applicable	Not applicable	- DC input power port

Electromagnetic emission Professional healthcare facility environment		
The CrystalView™ pump is designed for use in the electromagnetic environment specified below. The client or user of the CrystalView™ pump should ensure that it is used in such an environment.		
Phenomenon	Basic standard	Compliance level
Conducted and RF radiated voltage emissions	CISPR 11	Conducted emissions - Power supply 230Vac/50Hz (Phase and neutral between 150 kHz and 30 MHz) RF radiated voltage emissions (vertical and horizontal polarities; frequency between 30 MHz and 6 GHz) - Front side Group 1 Class B
Harmonic distortion	IEC 61000-3-2	AC power supply Class A
Voltage fluctuation and flicker	IEC 61000-3-3	AC power supply

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