DOUBLEFLO

Inflow/Outflow Pump

Instruction Manual

104611_H | April, 2021

REF.72205353

DAY TUBE SET

EN

TAGESSCHLAUCH-SET

DE

DAGSLANGENSET

NL

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IT

TUBULURE JOURNALIÈRE

FR

TUBULADURA DIARIA

ES

GÜNLÜK BORU SETİ

TR

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PL

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PT

DAGSSLANGESETT

NO

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SV


DAGSSLANGESÆT

DA

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FI

CE 0459

EN	DAY TUBE SET	4
DE	TAGESSCHLAUCH-SET	6
NL	DAGSLANGENSET	8
IT	SET DI TUBI GIORNALIERO	10
FR	TUBULURE JOURNALIÈRE	12
ES	TUBULADURA DIARIA	14
TR	GÜNLÜK BORU SETİ	16
PL	PRZEWÓD JEDNODNIOWY	18
PT	TUBAGEM DIÁRIA	20
NO	DAGSSLANGESETT	22
SV	DAGSLANGUPPSÄTTNING	24
DA	DAGSSLANGESÆT	26
FI	PÄIVÄLETKUSARJA	28

INDICATIONS FOR USE

This tubing is for daily use and delivers 0.9% saline solution (irrigation fluid) from a fluid bag to the Patient Tube Set. This device is used exclusively in conjunction with the Patient Tube Set (Ref. 72205354) and the DOUBLEFLO Pump (Ref. 72205352) during arthroscopic procedures.

INDICATION

Reserved for use by qualified health professionals. Must be used on 8-years-old or older pediatric and adults patients requiring an arthroscopic surgery, regardless of gender.

CONTRAINDICATION

The use of the device is prohibited when arthroscopy is warned against or disallowed. The known and possible side effects during arthroscopy procedure with a fluid management system are following :
Emphysema, hematoma, reflex sympathetic dystrophy, deep venous thrombosis (DVT) and fluid extravasations.

SUPPLY DETAILS

This tubing is sterilized with Ethylene Oxide. The tubing remains sterile as long as the package is unopened and undamaged before it gets to its final user. Re-usable for the surgical day when used with Ref. 72205354. For more details see warnings below and appropriate DOUBLEFLO Pump instruction manual (Ref. 104785).

WARNING

Do not press on the "START/STOP" button of the DOUBLEFLO Pump before setting up all tubing.

WARNING

The Day Tube Set can be re-used for multiple procedures during the same surgical day only when used in conjunction with the Patient Tube Set Ref. 72205354. The Patient Tube Set must be changed for each patient.

WARNING

1. Do not use if the package is damaged or opened.
2. Sterility is preserved only if proper setup and tube changing procedures have been followed according to the DOUBLEFLO Pump instruction manual (Ref. 104785).
3. Refer to and read the DOUBLEFLO Pump instruction manual (Ref. 104785) for proper use.
4. Re-sterilization: This product has not been designed to be re-sterilized. Reprocessing may lead to changes in material characteristics such as deformation and material degradation which may impact the strength of the device and compromise device performance. These risks may potentially affect patient safety.
5. Only use with a Patient Tube Set with check valve Ref. 72205354.
6. Replace the tubing in case of malfunction or variation in performance.
7. Be careful when inserting or removing the blue connector: do not insert nor remove the blue connector with pressure inside to avoid to deform the white silicone membranes of the pressure sensors of the pump.

Any serious incidents that have occurred in connection with these devices must be notified to the manufacturer and the competent authority of the Member State in which the user is based.

INSTRUCTIONS FOR USE (to be performed by the circulating nurse)

1. Open the package using aseptic technique.

The pump should be on standby when setting up this tubing.

2. Open the outer package by using the opening tab
3. Remove the tube set from the inner package and close the two clamps
4. Remove all the caps on the spikes and connect the spikes to the bags of saline solution

Insert the blue connector only if the membranes are not pressurized. Maintain the blue connector vertically throughout insertion.

Do not touch the white silicone membranes of the pressure sensors of the pump.

6. Insert the blue connector from the tubing into the bracket with blue indicator (right opening)
7. Place and center the tube around the irrigation pump head

Do not twist the tube.

8. Pull the tube tight through the left opening of the black bracket for the irrigation side (left opening)
9. Insert the white connector from the tubing into the black bracket of the irrigation side (left opening)
10. Connect the blue end ("Luer-Lock") to the Patient Tube Set (Ref. 72205354); refer to the Patient Tube Set instruction manual (Ref. 104612)

ONCE ALL TUBE SETS ARE INSTALLED

- I. The circulating nurse opens one clamp of the Day Tube Set (Ref. 72205353)
- II. The scrub nurse opens the arthroscopy sheath / cannula stopcock
- III. The circulating nurse presses on the "START/STOP" button of the DOUBLEFLO Pump (Ref. 72205352) when ready to begin to fill the irrigation tubing
- IV. When the irrigation tubing has been purged of air, the circulating presses again on the "START/STOP" button of the DOUBLEFLO pump (Ref. 72205352)

END OF OPERATING DAY

- I. Switch off the fluid management and ensure that the membranes of the blue connector are not pressurized
- II. Close all clamps, remove the spikes from the saline solution bags, then open the clamp near the red "Luer-Lock" and disconnect the white and blue connector from the brackets

For more details, refer to the DOUBLEFLO pump instruction manual (Ref. 104785).

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SMITH & NEPHEW, Inc.
HEMODIA SAS