DOUBLEFLO

Inflow/Outflow Pump

Instruction Manual

104612_H | April, 2021

REF.72205354

PATIENT TUBE SET

EN

PATIENTENSCHLAUCH-SET

DE

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NL

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IT

TUBULURE PATIENT

FR

TUBULADURA PACIENTE

ES

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TR

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PL

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PT

SETT MED PASIENT SLANGE

NO

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SV

PATIENTSLANGESÆT

DA

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FI















CE 0459

To Day Tube Set

FLOW

To Arthroscopy Sheath

PATIENT TUBE SET

REF. 72205354

	To Day Tube Set	FLOW	To Arthroscopy Sheath
EN	PATIENT TUBE SET		6
DE	PATIENTENSCHLAUCH-SET		8
NL	PATIËNTSLANGENSET		10
IT	SET DI TUBI PAZIENTE		12
FR	TUBULURE PATIENT		14
ES	TUBULADURA PACIENTE		16
TR	HASTA BORU SETİ		18
PL	PRZEWÓD PACJENTA		20
PT	TUBAGEM DO PACIENTE		22
NO	SETT MED PASIENT SLANGE		24
SV	PATIENTSLANGUPPSÄTTNING		26
DA	PATIENTSLANGESÆT		28
FI	POTILASLETKUSARJA		30

INDICATION FOR USE

The device is intended to be used exclusively in conjunction with the DOUBLEFLO Pump (Ref. 72205352) and delivers 0.9% saline solution (irrigation fluid) from the Day Tube Set (Ref. 72205353) to the arthroscopy sheath / cannula during arthroscopic procedures.

INDICATION

Reserved for use by qualified health professionals. Must be used on 8-years-old or older pediatric and adults patients requiring an arthroscopic surgery, regardless of gender.

CONTRAINDICATION

The use of the device is prohibited when arthroscopy is warned against or disallowed. The known and possible side effects during arthroscopy procedure with a fluid management system are following :

Emphysema, hematoma, reflex sympathetic dystrophy, deep venous thrombosis (DVT) and fluid extravasations.

SUPPLY DETAILS

Single use device, sterilized with Ethylene Oxide. The tubing remains sterile as long as the package is unopened and undamaged before it gets to its final user.

WARNING

Do not press on the "START/STOP" button of the DOUBLEFLO Pump before setting up all tubing.

WARNING

To properly connect the male Luer connector on the female Luer connector of the arthroscopy sheath / cannula, the ring of the male Luer connector is moved backward (Fig. VI.a and VI.b); then the male cone of the male Luer connector is introduced into the female cone of the female Luer connector to the stop and do a quarter turn (Fig. VII.a and VII.b); then the ring is moved forward and screwed onto the screw thread of the female Luer connector (Fig. VIII.a and VIII.b).

WARNING

1. Do not use if the package is damaged or opened.
2. Sterility is preserved only if proper setup and tube changing procedures have been followed according to the DOUBLEFLO Pump instruction manual (Ref. 104785).
3. Refer to and read the DOUBLEFLO Pump instruction manual (Ref. 104785) for proper use.
4. This product is for single-use only. It has not been designed to be re-used/re-sterilized. Reprocessing may lead to changes in material characteristics such as deformation and material degradation which may impact the strength of the device and compromise device performance. Reprocessing of single use devices can also cause cross-contamination leading to patient infection. These risks may potentially affect patient safety.
5. Never attach one check valve to another. This could cause loss of tubing performance or a risk of cross-contamination.
6. Replace the tubing in case of malfunction or variation in performance.

Any serious incidents that have occurred in connection with these devices must be notified to the manufacturer and the competent authority of the Member State in which the user is based.

INSTRUCTIONS FOR USE

1. Open the package using aseptic technique

The pump should be on standby when setting up the tubing.

2. The circulating nurse opens the outer package by using the opening tab
3. The scrub nurse removes the tube set from the inner package and ensures that the clamp is opened (otherwise, open the clamp)
4. The scrub nurse ensures that all “Luer-Locks” are tight, then hands the tubing end with the blue cap to the circulating nurse
5. The circulating nurse removes the blue cap from the end of the Patient Tube Set (Fig. I. D) and the Day Tube Set (Ref. 72205353) (Fig. I. C)
6. The circulating nurse connects the Patient Tube Set (Fig. II. D) to the Day Tube Set (Ref. 72205353) (Fig II. C) and tightens the ring to lock both Luer connectors (Refer to Fig VI.a, VII.a, VIII.a and warnings to properly connect the male Luer connector)
7. The scrub nurse connects the red “Luer-Lock” from the Patient Tube Set to the arthroscopy sheath / cannula (Refer to Fig. VI.b, VII.b, VIII.b and warnings to properly connect the male Luer connector)

ONCE ALL TUBE SETS ARE INSTALLED

- I. The circulating nurse opens one clamp of the Day Tube Set (Ref. 72205353)
- II. The scrub nurse opens the arthroscopy sheath / cannula stopcock
- III. The circulating nurse presses on the “START/STOP” button of the DOUBLEFLO Pump (Ref. 72205352) when ready to begin surgery to fill the irrigation tubing
- IV. When the irrigation tubing has been purged of air, press again on the “START/STOP” button of the DOUBLEFLO pump (Ref. 72205352)

END OF THE OPERATION

Flush the irrigation system after each patient procedure prior to disconnection.

- I. The circulating nurse closes the opened clamp on the Day Tube Set (Ref. 72205353) prior to disconnection
- II. The scrub nurse closes the clamp on the Patient Tube Set
- III. The scrub nurse disconnects the red “Luer Lock” from the Patient Tube Set to the arthroscopy sheath / cannula
- IV. The circulating nurse disconnects the blue “Luer Lock” (Fig. III. E) under the check valve and dispose of this tubing in a dedicated site (part of the Patient Tube Set)

The check valve will remain connected to the Day Tube Set (Ref. 72205353) until the next arthroscopic case.

START OF THE NEXT OPERATION

- I. Use aseptic technique
- II. The circulating nurse disconnects the blue “Luer Lock” to the Day Tube Set (Ref. 72205353) (Fig. IV. C/D) and dispose of the tube with check valve in a dedicated site

Must take off her gloves and wear new ones before connecting the new Patient Tube Set.

- III. Then the circulating nurse immediately connects the blue Luer-lock (Fig. V. D) of the new Patient Tube Set to the blue connector (Fig. V. C) of the Day Tube Set (Ref. 72205353) as per step 6 in the Instructions for use (Refer to Fig. VI.a, VIIa, VIIIa and warnings to properly connect the male Luer connector)
- IV. Refer to the step 7 of “INSTRUCTION FOR USE” to IV of “ONCE ALL TUBE SETS ARE INSTALLED” above

END OF OPERATING DAY

- I. Switch off the fluid management and ensure that the membranes of the blue connector are not pressurized
- II. Close all clamps, remove the spikes from the saline solution bags, then open the clamp near the red “Luer-Lock” and disconnect the white and blue connector from the brackets

For more details, refer to the DOUBLEFLO pump instruction manual (Ref. 104785).

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SMITH & NEPHEW, Inc.
HEMODIA SAS