DOUBLEFLO

Inflow/Outflow Pump

Instruction Manual

104613_H | April, 2021

REF.72205355

INFLOW TUBE SET

EN

IRRIGATIONSSCHLAUCH-SET

DE

IRRIGATIESLANGENSET

NL

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IT

TUBULURE D'IRRIGATION

FR

TUBULADURA DE IRRIGACIÓN

ES

SULAMA BORU SETİ

TR

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PL

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PT

SETT MED TILSTRØMMINGSSSLANGE

NO

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SV

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




DA

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
FI

CE 0459

INFLOW TUBE SET

REF. 72205355

EN	INFLOW TUBE SET	6
DE	IRRIGATIONSSCHLAUCH-SET	8
NL	IRRIGATIESLANGENSET	10
IT	SET DI TUBI D'IRRIGAZIONE	12
FR	TUBULURE D'IRRIGATION	14
ES	TUBULADURA DE IRRIGACIÓN	16
TR	SULAMA BORU SETİ	18
PL	PRZEWÓD NAWADNIAJĄCY	20
PT	TUBAGEM DE IRRIGAÇÃO	22
NO	SETT MED TILSTRØMMINGSSLANGE	24
SV	INFLÖDESSLANGUPPSÄTTNING	26
DA	INDLØBSSLANGESÆT	28
FI	TULOVIKTAUSLETKUSARJA	30

INDICATIONS FOR USE

The device delivers exclusively 0.9% saline solution (irrigation fluid) from the bags to the arthroscopy sheath / cannula. This device is used exclusively in conjunction with the DOUBLEFLO Pump (Ref. 72205352) during arthroscopic procedures.

INDICATION

Reserved for use by qualified health professionals. Must be used on 8-years-old or older pediatric and adults patients requiring an arthroscopic surgery, regardless of gender.

CONTRAINDICATION

The use of the device is prohibited when arthroscopy is warned against or disallowed. The known and possible side effects during arthroscopy procedure with a fluid management system are following :

Emphysema, hematoma, reflex sympathetic dystrophy, deep venous thrombosis (DVT) and fluid extravasations.

SUPPLY DETAILS

Single use device, sterilized with Ethylene Oxide. The tubing remains sterile as long as the package is unopened and undamaged before it gets to its final user.

WARNING

Do not press on the "START/STOP" button of the DOUBLEFLO Pump before setting up all tubing.

WARNING

To properly connect the male Luer connector on the female Luer connector, the ring of the male Luer connector is moved backward (Fig. I.a and I.b); then the male cone of the male Luer connector is introduced into the female cone of the female Luer connector to the stop (Fig. II.a and II.b); then the ring is moved forward and screwed onto the screw thread of the female Luer connector (Fig. III.a and III.b).

WARNING

1. Do not use if the package is damaged or opened.
2. Sterility is preserved only if proper setup and tube changing procedures have been followed according to the DOUBLEFLO Pump instruction manual (Ref. 104785).
3. Refer to and read the DOUBLEFLO Pump instruction manual (Ref. 104785) for proper use.
4. This product is for single-use only. It has not been designed to be re-used/re-sterilized. Reprocessing may lead to changes in material characteristics such as deformation and material degradation which may impact the strength of the device and compromise device performance. Reprocessing of single use devices can also cause cross-contamination leading to patient infection. These risks may potentially affect patient safety.
5. Replace the tubing in case of malfunction or variation in performance.
6. Be careful when inserting or removing the blue connector: do not insert nor remove the blue connector with pressure inside to avoid to deform the white silicone membranes of the pressure sensors of the pump.

Any serious incidents that have occurred in connection with these devices must be notified to the manufacturer and the competent authority of the Member State in which the user is based.

INSTRUCTIONS FOR USE

1. Open the package using aseptic technique.

The pump should be on standby when setting up this tubing.

2. The circulating nurse opens the outer package by using the opening tab
3. The scrub nurse removes the tube set from the inner package and ensures that the clamp close to the tube dedicated to the arthroscopy sheath / cannula is opened (otherwise, open the clamp)
4. The scrub nurse keep the red end ("Luer-Lock") and hands the rest of the tubing with the two spikes to the circulating nurse
5. The circulating nurse closes the two clamps on the tubes dedicated to the connection with the saline bags, removes all the caps on the spikes and connect the spikes to the bags of saline solution

Insert the blue connector only if the membranes are not pressurized.
Do not touch the white silicone membranes of the pressure sensors of the pump.

Avoid to touch the pressure sensors

6. The circulating nurse inserts the blue connector of the tubing into the bracket with blue indicator (right opening)
7. The circulating nurse places and centers the tube around the roller head.

Do not twist the tube.

8. The circulating nurse pulls the tube tight through the left opening of the black bracket for the irrigation side (left opening)
9. The circulating nurse inserts the white connector from the tubing into the black bracket of the irrigation side (left opening)
10. The scrub nurse connects the red end ("Luer-Lock") of the Inflow Tube Set to the arthroscopy sheath (Refer to Fig. I.a, II.a, III.a and warnings to properly connect the male Luer connector)

ONCE ALL TUBE SETS ARE INSTALLED

- I. The scrub nurse opens the arthroscopy sheath / cannula stopcock
- II. The circulating nurse opens one of the clamps on the proximal end of the Inflow Tube Set (Ref. 72205355) connected to the saline bags.
- III. The circulating nurse presses "START/STOP" button of the DOUBLEFLO Pump (Ref. 72205352) when ready to begin to fill the irrigation tubing.
- IV. When the irrigation tubing has been purged of air, the circulating nurse presses again on the "START/STOP" button of the DOUBLEFLO pump (Ref. 72205352)

END OF OPERATING DAY

- I. Switch off the fluid management and ensure that the membranes of the blue connector are not pressurized
- II. Close all clamps, remove the spikes from the saline solution bags, then open the clamp near the red "Luer-Lock" and disconnect the white and blue connectors from the brackets

For more details, refer to the DOUBLEFLO pump instruction manual (Ref. 104785).

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SMITH & NEPHEW, Inc.
HEMODIA SAS