DOUBLEFLO

Inflow/Outflow Pump

Instruction Manual

104614_H | April, 2021

REF.72205356

OUTFLOW TUBE SET

ASPIRATIONSSCHLAUCH-SET

ASPIRATIESLANGENSET

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TUBULURE D'ASPIRATION

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EMME BORU SETİ

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EN

DE

NL

IT

FR

ES

TR

PL

PT

NO

SV


DA

FI

CE 0459

OUTFLOW TUBE SET

REF. 72205356

EN	OUTFLOW TUBE SET	4
DE	ASPIRATIONSSCHLAUCH-SET	6
NL	ASPIRATIESLANGENSET	8
IT	SET DI TUBI DI ASPIRAZIONE	10
FR	TUBULURE D'ASPIRATION	12
ES	TUBULADURA DE ASPIRACIÓN	14
TR	EMME BORU SETİ	16
PL	PRZEWÓD SSĄCY	18
PT	TUBAGEM DE ASPIRAÇÃO	20
NO	SETT MED UTSTRØMMINGSSSLANGE	22
SV	UTFLÖDESSLANGUPPSÄTTNING	24
DA	UDLØBSSLANGESÆT	26
FI	LÄHTÖVIRTAUSLETKUSARJA	28

INDICATIONS FOR USE

The device is intended to be used exclusively in conjunction with the DOUBLEFLO Pump (Ref. 72205352) during arthroscopic procedure and is used to remove the fluid in the joint via the cannula or a shaver.

INDICATION

Reserved for use by qualified health professionals. Must be used on 8-years-old or older pediatric and adults patients requiring an arthroscopic surgery, regardless of gender.

CONTRAINDICATION

The use of the device is prohibited when arthroscopy is warned against or disallowed. The known and possible side effects during arthroscopy procedure with a fluid management system are following :

Emphysema, hematoma, reflex sympathetic dystrophy, deep venous thrombosis (DVT) and fluid extravasations.

SUPPLY DETAILS

Single use device, sterilized with Ethylene Oxide. The tubing remains sterile as long as the package is unopened and undamaged before it gets to its final user.

WARNING

Do not press on the "START/STOP" button of the DOUBLEFLO Pump before setting up all tubing.

WARNING

1. Do not use if the package is damaged or opened.
2. Sterility is preserved only if proper setup and tube changing procedures have been followed according to the DOUBLEFLO Pump instruction manual (Ref. 104785).
3. Refer to and read the DOUBLEFLO Pump instruction manual (Ref. 104785) for proper use.
4. This product is for single-use only. It has not been designed to be re-used/re-sterilized. Reprocessing may lead to changes in material characteristics such as deformation and material degradation which may impact the strength of the device and compromise device performance. Reprocessing of single use devices can also cause cross-contamination leading to patient infection. These risks may potentially affect patient safety.
5. Replace the tubing in case of malfunction or variation in performance.

Any serious incidents that have occurred in connection with these devices must be notified to the manufacturer and the competent authority of the Member State in which the user is based.

INSTRUCTIONS FOR USE

1. Open the package using aseptic technique.

The pump should be on standby when setting up this tubing.

2. The circulating nurse opens the outer package by using the opening tab
3. The scrub nurse removes the tube set from the inner package and ensures that the clamp is opened (otherwise, open the clamp)
4. The scrub nurse hands the tube end marked "WASTE" to the circulating nurse
5. The circulating nurse inserts the "WASTE" end into the waste collection system
6. The circulating nurse inserts the orange connector of the tubing into the bracket with orange indicator (left opening)
7. The circulating nurse places and centers the tubing around the suction pump head.
8. The circulating nurse inserts the white connector from the tubing into the black bracket (right opening)

Do not twist the tube.

9. The circulating nurse inserts the "OUTFLOW CANNULA" tube into the right side of the pinch valve
10. The circulating nurse inserts the "SHAVER" tube into the left side of the pinch valve

Make sure the DOUBLEFLO pump is positioned such that the "OUTFLOW CANNULA" tube and "SHAVER" tube hang vertically and the shelf does not interfere.

11. The scrub nurse connects the tubing end marked "OUTFLOW CANNULA" onto the cannula
12. The scrub nurse connects the tubing end marked "SHAVER" onto the shaver device

ONCE ALL TUBE SETS ARE INSTALLED

- I. The circulating nurse opens one of the clamps of the Day Tube Set (Ref. 72205353) or Inflow Tube Set connected to the saline bags.
- II. The circulating nurse presses the "START/STOP" button of the DOUBLEFLO Pump (Ref. 72205352) when ready to begin to fill the irrigation tubing.
- III. The scrub nurse opens the arthroscopy sheath / cannula stopcock.

END OF OPERATING DAY

- I. Drain the tubing and disconnect the accessories
- II. Disconnect the white connector from the bracket, the "OUTFLOW CANNULA" and "SHAVER" end marked from the pinch valve then the orange connector from the bracket

For more details, refer to the DOUBLEFLO pump instruction manual (Ref. 104785).

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SMITH & NEPHEW, Inc.
HEMODIA SAS