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Instructions for use - Dialysis Connectors

Read instructions carefully before use

EN

Packaging and dimensions

Flexible PA (polyamide)/PE (polyethylene) and PET (polyethylene terephthalate)/PE (polyethylene) blister pack containing one or two dialysis connectors of variable length. The number of connectors per blister pack (one or two lines) and the length of the tubing may vary depending on the product.

Composition

The dialysis connectors consist of tubing comprising:

- Male Luer lock connector with a rotating collar **A** and female LL cap **B**
- Tubes **C**
- Female LL connector **D** + male LL cap **E**
- Clamps **F**

Indication and clinical use

Purpose of the device: the connectors are used during haemodialysis. They are designed to extend blood lines. They are connected to the venous access device (catheter or needle) to transport blood during the entire dialysis session (approximately 4 hours).


Compatible devices: simple dialysis connectors are designed to extend blood lines during haemodialysis performed using a fistula (needle) or catheter.

Instructions for connecting Luer Lock fittings

Haemodialysis is carried out according to current protocols

Clamps are used to close and open lines according to protocol requirements
Colour:

- Red is for arterial lines
- Blue is for venous lines
- White can be used for both arterial and venous lines


Step 1:

Slide the collar of the male Luer lock fitting backwards and HOLD it in this position


Step 2:

Holding the collar back, insert the male tip of the male Luer lock fitting into the female tip of the female Luer fitting until it will not go further.


Step 3:

Rotate a quarter turn, still holding the collar back and ensuring the male and female tips are fully connected. The tip-to-tip connection is now complete and the fittings are mechanically linked.

Step 4:

Slide the male Luer lock fitting collar forwards and then screw it onto the female Luer lock fitting until it will not turn any more.

Warnings and precautions for use

- Please read all the information in these Instructions for Use carefully before using the dialysis connectors.
- The tubing must only be used for its intended purpose. Healthcare professionals are personally responsible for using these dialysis connectors in accordance with their intended use.
- The device is sterilised by gamma irradiation.
- Check the integrity of the packaging before use.
- Do not use the device if the packaging is damaged or open, as the sterility of the device may be compromised.
- The device is not designed to be re-sterilised. Re-sterilisation could alter the characteristics of the materials, leading to their deformation and degradation, which may affect the rigidity of the device and compromise its performance.
- Single-use device. Reprocessing and reuse of the single-use device may result in infection or cause other harm to the patient.
- Check that the connectors are free from kinks during fitting and/or use of the device.
- Check that the connectors are tightly connected before starting dialysis.
- Check the expiry date on the packaging or product label. Do not use if the product has expired.
- HEMODIA expressly prohibits modifications of any kind to its connectors. Any modification or deviation from the instructions provided herein shall invalidate the liability of HEMODIA.

As a medical device manufacturer, HEMODIA cannot be held responsible for any incident resulting from improper use of its line.

Contraindications

- Do not use in children under 18 years of age
- Do not re-sterilise
- Do not reuse

Conditions for storage and disposal

Storage: Store in a dry place.
Store protected from light.

Disposal: refer to hospital waste disposal methods

Medical device vigilance

Any serious incident involving this device must be reported to the manufacturer or to the competent health authority. By reporting serious incidents, you are helping to provide more information about the safety of the device.
<https://ansm.sante.fr/>

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- Publication date: September 2023
- Code: 105332_A

Label Symbols Key

	CE 0459 This device is labelled with the CE Marking, which was obtained from the notified body GMED in 2004
	Hemodia Manufacturer's Reference
	Batch number
	Expiry date
	Latex free
	Note: please refer to the instructions for use for important safety information, such as warnings and precautions.
	Manufacturer
	Store in a dry place away from moisture
	Keep away from direct sunlight
	Sterilised by Gamma irradiation
	Do not reuse
	Check the integrity of the packaging before use
	Single sterile barrier system
	Refer to instructions for use
	Quantity per box

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Hemodia SAS