

**Praxiject™ SF Saline Normale**  
**Prefilled Syringe with 0.9% Sodium Chloride**

**INDICATIONS AND CLINICAL USE**

For vascular access devices flushing only. May be placed on a sterile field. For use by healthcare professionals trained in lock / flush procedures.

**CONTRAINDICATIONS**

Not for dry product reconstitution or for medication dilution.

**WARNINGS AND PRECAUTIONS**

Use aseptic technique. Single use only. Do not reuse. Risk of contamination if reused. Not for direct intravenous injection. Do not use if package is opened or damaged. Do not use if the cap on the syringe is not intact. Visually inspect the contents of each prefilled syringe for clarity, particulate matter, precipitate, discoloration, and leakage prior to use. Do not use if any of the aforementioned is observed. Do not allow air to be trapped in fluid path. Do not resterilise. Sterile unless package opened or damaged. Serious incidents that occur in relation to the device should be reported to the manufacturer and the competent authority of the user/patient's Member State.

**ADVERSE REACTIONS**

No known adverse reactions when the product is used as indicated.

**INSTRUCTIONS FOR USE**

**USE ASEPTIC TECHNIQUE.**

1 Aseptically open package containing the 0.9% Sodium Chloride prefilled syringe. Visually inspect product for particulate matter and discoloration prior to use according to WARNINGS AND PRECAUTIONS.

2 Remove protector as required. Expel air from syringe. Do not allow air to be trapped in fluid path. Use only with compatible luer lock connectors.

3 In accordance with the institution's protocol, flush the lumen of the vascular access device (VAD) with 0.9% Sodium Chloride solution to remove any blood, medication or other substance remaining in the VAD.

4 Remove and discard unused portions and empty syringes according to the institution's biohazardous waste policy.

**ACTION**

0.9% Sodium Chloride Injection USP solution is a sterile, aqueous solution having approximately the same osmotic pressure and composition as extracellular fluids. It is non-irritating to tissues. It is used to flush vascular access devices in order to maintain catheter patency and to prevent contact between incompatible medications or fluids.

**PRODUCT DESCRIPTION AND PACKAGING**

Each syringe is supplied as a clear, latex-free, plastic prefilled syringe, packaged in a heat-sealed pouch. Each syringe contains 0.9% Sodium Chloride Injection USP with no preservatives.

REF	Description	Packaging
1100-4	10mL 0.9% Sodium Chloride Injection USP in a 10mL syringe	One (1) 10mL syringe packaged in a heat sealed pouch, 100 pouches <a href="http://www.merckusj.com">www.merckusj.com</a>

**STORAGE AND STABILITY**

0.9% Sodium Chloride Injection USP should be stored between 15°C and 30°C (59°F and 86°F), protected from direct sunlight and freezing. Product may be stored as indicated by the expiration date on the product package.

Doc. No.: MBX3705x-4 (EU)\_v1

REV Date: 2023-01-20