

INFORMATION ABOUT ASSEMBLIES

USAGE: This assembly is single-use and intended for a single patient.

Warnings and Precautions for use:

- Check the integrity of the packaging before use, do not use the device if the packaging has been damaged or opened before use.
- Do not resterilize.
- Do not reuse : Reuse of the assembly may cause infection or harm to the patient
- Check the expiry date on the product package or label. Do not use if the product has expired.
- PLEASE NOTE that your pack may contain one or more items that require precautionary measures specified on the product or accompanying documentation.

STORING CONDITIONS : Please refer to the product label.

Medical devices vigilance:

Any serious incident occurring in connection with this assembly must be notified to the manufacturer or to the relevant health authority. By reporting serious incidents, you help to provide more information on the safety of the device.

Please consult the special instructions for use on the back.

Creation date : MARCH 15, 2024

Update Rev.1 : APRIL 30th, 2024

Page 1/2



ASSEMBLED BY :
Hemodia SAS
85 rue du Chêne vert
31670 LABEGE- FR

Information AA21

Component	Information
HAND TOWEL	For cleaning of equipment or hand hygiene.
WASHCLOTH	Product for personal hygiene and patient comfort.
WASTE BAG	For non-sharp, non-prickly or unsoiled waste
SYRINGE	The syringes in the set are made of polypropylene and can contain polycarbonate, silicone, isoprene and elastomer. For any use with drug substances, please refer to the usage notice of the substance.
TRAYS	Trays are packaging elements, which are not intended to contain injectable solutions.
GALLIPOTS	The cups are made of polypropylene. They are intended only for containing an antiseptic liquid or physiological serum. They are not intended for containing implantable devices or injectable preparations. The graduations do not have a measuring function; they only give an indication of volume.
SAVONGE	The savonge is intended for the personal hygiene of the patients. Not recommended for use in children under two years of age.
MASKS	The medical face masks comply with standard EN 14683:2019. Unless otherwise stated on the label, the masks must be Type II (bacterial filtration efficiency ≥98%). Wash hands. Secure ear loops behind the ears, with the blue side on the outside. Make sure it covers the nose and mouth. Do not pull the mask down or up. Avoid touching the mask.

Creation date : MARCH 15, 2024

Update Rev.1 : APRIL 30th, 2024

Page 2/2