

DISSECTING FORCEPS

Instruction for use AA76

INTENDED USE

Used to hold tissue in place (wound edges) when applying sutures, to move tissue out of the way during surgery or to move or remove dressings without using the hands or fingers

TESTING AND INSPECTION:

Inspect for rust, pits, cracks, pores, and dents before use.

WARNINGS:

Don't use for second time

Don't use if you observe sharp edges

Product material is conductive : take away from electric source.

Do not use the wrong instrument, as it may cause problems during operation. (consigne de sécurité)

STORAGE & TRANSPORT CONDITIONS :

Ambient temperature range : 0°C to 50°C (32°F to 122°F)- Relative humidity range : 30% to 75% - Atmospheric pressure range : 700 hPA to 1060 hPA

RECYCLE:

When processing medical devices please follow your local government Health & Safety procedures.

 SURGIKARE
Toor Abad Daska Road,
Sialkot, Punjab - PK

 Obelis
Bd Général Wahis, 53
B-1030 Brussels - BE

 **€1639** 