


## ONE USE HOLDER SAFETY INFORMATION



THIS HOLDER IS FOR SINGLE USE ONLY. USE ONCE AND DISCARD. MULTIPLE USES OF THE NEEDLE HOLDER MAY RESULT IN DAMAGED THREADS AND CAUSE THE NEEDLE TO UNTHREAD OR SPIN OUT.

Handle all biologic samples and blood collection needles in accordance with the policies and procedures of your facility. Obtain appropriate medical attention in the event of any exposure to biologic samples (for example, through a puncture injury) since samples may transmit viral hepatitis, HIV (AIDS), or other infectious diseases. Intended for use by trained healthcare professionals only. EU Only: Users should report any serious incident related to the device to the Manufacturer and National Competent Authority. Outside EU: Contact your local BD representative for any incident or inquiry related to this device.

Technical Services and Support: In the United States, contact BD at 1-800-638-8663 OPT 3 or [bd.com](http://bd.com).

For regions outside of the United States, contact your local BD representative or [bd.com](http://bd.com).

For U.S. patents that may apply, see [bd.com/patents](http://bd.com/patents).

2023-08 VDM60736(01) 500072214

## Instruction for use AA98\_EN

 Becton, Dickinson and Company,  
1 Becton Drive, Franklin Lakes,  
New Jersey 07417 - 1885 - USA



Becton Dickinson Ireland Ltd,  
Donore Road, Drogheda  
Co. Louth, A92 YW26 - IE



Becton Dickinson Distribution  
Center NV, Laagstraat 57,  
9140 Temse, Belgium

Internal drafting date 13/01/2026