

AQUA VISION[®]

Fluid management pump
for arthroscopic procedures



Shoulder



Knee




Small joint



Hip






AQUA VISION®

The fluid management pump with integrated shaver interface

The AQUA VISION® is intended to be used during arthroscopic procedures by orthopedic surgeons and operating room nurses


- A pump and its accessories developed in collaboration with surgeons
- Simple, user-friendly operation that allows for rapid mastery of the device
- Integrated shaver interface, which allows the use of multiple types of shavers on the market



The Hemodia experience

A recognised knowledge

- Hemodia has developed leading-edge expertise in fluid management, both for pumps and tubing
- Its recognised know-how provides solutions for equipment, consumables and pump maintenance
- All AQUA VISION® pumps are assembled in our workshop in Toulouse (France)



Adopt our products and services



Quality

Hemodia manufactures with a constant concern for improving its products and in compliance with standards



Reliability

Hemodia has a fully in-house process and several integrated production sites



Customer service

The customer service and the Hemodia sales teams are at your side in the field to meet your expectations

Tubing adapted to the pump

Our products



Daily irrigation tubing

The daily use device allows irrigation fluid to be delivered from the bag to the intermediate irrigation tubing



Intermediate irrigation tubing

This device is designed to deliver the irrigation fluid from the daily irrigation tubing to the arthroscope



Suction tubing

Used for fluid extraction from the joint via cannula or shaver

Tubing set up on the AQUA VISION® pump

The pump is designed to perfectly meet the demands of IBODEs, thanks to which the set up is:




Simple



Intuitive



Secure



Interfaces

- **The Integrated Shaver Interface (ISI), the shaverinterfaceclip** works with all shaver brands. The pump senses the operation of the handpiece and automatically adjusts the shaver/cannula suction
- **The Foot Control Interface** works with Smith & Nephew®, ConMed®, Arthrex®, Eberle® and Karl Storz®



About Hemodia

Founded in 1985, Hemodia has developed expertise in medical devices. Based on core business, medical tubing and polymer components, Hemodia has rapidly become a major player on the custom-made care kits market. From dialysis to anesthesia sets, Hemodia adjusts its offer to every request and speciality : oncology, infusion, hygiene...

Thanks to the concept of internalisation of conception, production, sterilisation, storage and distribution – Hemodia anticipates the customer need and answers all requirements concerning deliveries once orders have been placed. With production in France and Tunisia it ensures total control on supplies and equipment.

Thanks to Quality Management System, Hemodia stands by your side by reducing the risk of infections, improving compliance and best practice, avoiding wastage and creating value.

The high level of requirements guarantees both professionals and patients the best quality products.

Trusting in Hemodia is choosing knowledge and recognised expertise.

Technical specifications

Equipment

Reference	Description
ARTH-PUMP-A	AQUA VISION® arthroscopy pump
ARTH-PED2-A	2-way pedal
ARTH-PED4-A	4- way pedal
ARTH-ISHAU	Shaver interface clip
ARTH-FCI01	Foot Control Interface for Smith & Nephew / Dyonics shaver consoles
ARTH-FCI02	Foot Control Interface for ConMed / Linvatec shaver consoles
ARTH-FCI03	Foot Control Interface for Arthrex shaver consoles
ARTH-FCI05	Foot Control Interface for Karl Storz shaver consoles
103773	AC power cord for Europe
103822	AC power cord for United Kingdom
103899	AC power cord for Switzerland
104270	AC power cord for Canada
103777	Black carrying case

Consumables

Reference	Description	Units per box
AV4503	Day set	30
AV4509	Patient set with check valve	30
AV1102	Intermediary tubing with check valve	50



For more informations please contact :

Customer support

+33 (0)5 61 00 71 81
customer.support@hemodia.com



Hemodia SAS France

85 rue du Chêne vert
31670 Labège - France

Your local representative



Class IIa consumables and pump - Other Class I equipment (except carrying case).
Please read the instructions for use accompanying these medical devices carefully before use.

