

Catalogue

Custom-made care kits



Table of contents



General information

p.2



Catalogue

p.3



Anesthesia kits

- Spinal
- Epidural
- Nerve blocks

p.5



Dialysis kits

- CAPD
- Catheter
- Fistula ON/OFF kits

p.9



Care kits

- Dressing
- Suturing
- Sponge stick

p.17



Infusion kits

- Central venous access CVC
- Peripheral venous catheter PVC
- Peripheral venous neuraxial

p.29



Central venous access fitting kits

- Implantable chamber
- Catheter
- PICCLine
- MIDLine

p.39



Catheterisation kits

- Drain
- Tubes
- Urinary catheterisation

p.43

General information



Hemodia

Thanks to its annual production of 20 millions units, Hemodia is the French leader in the manufacturing of customised single-use care kits.

Customised compositions adapted to your needs

For more than 30 years, Hemodia has developed its industrial tools and know-how.

These customised care kits are class I to IIa medical devices, CE marked by Hemodia SAS. Their class is determined according to their composition. For class Is and IIa, the notified body GMED (CE0459) has been involved.

Read carefully the instructions in the package insert and on the label.

With a wide range of components, our care kits meet your specifications and current regulatory requirements.

Hemodia's sales teams are at your service and support you in the development of technical and logistics solutions.





Customised care kits

On the following pages, you will find an overview of our range of products.

For each indicator, a component list is given as a guideline to help you build your product.


Our packaging solutions allow many different configurations to best adapt to your request.

For any not listed components, a feasibility study can be carried out.

For some components, Hemodia is not the legal manufacturer. View list P46 and P47.



For all non catalogue related requests, please contact our sales team.



Anesthesia kits

- Spinal
- Epidural
- Nerve blocks

Loco-Regional anesthesia kit

Protective and hygiene components :

- Gown or surgical gown
- Mob cap
- Face mask
- *Surgical brush*
- Hand towel
- Table drape
- Probe cover
- Patient drape with aperture and adhesives

Patient prepping components :

- Sponge stick
- Swabs
- Gallipots
- Forceps
- Sponge forceps

Technical components** :

- Blunt fill needle
- Hypodermic needle
- Spinal needle*
- Caps
- Extension line
- IV catheter
- Syringes
- *LOR syringes*
- Filters

Others:

- Scissors
- Tourniquet
- Ultrasound transmission gel
- Dressings
- Clamping forceps
- Adhesive strips
- DASRI bag

Italic text refers to components currently under sourcing (not yet available)

*These products are class III medical devices, CE marked by BEXTON DICKINSON S.A.; the notified body Agencia española de medicamentos y productos sanitarios (CE0318) has intervened.

**Following ISO 80369 a new standard connection is used for regional anaesthesia: NRFit, some NRFit components are already available (syringes and fill needles) please contact customer support for more information



This composition is an example, we adapt to your protocol.



Dialysis kits

- CAPD
- catheter
- Fistula

Dialysis Fistula ON kit

Protective and hygiene components :

- Gloves
- Gowns
- Face mask
- Table drape
- Patient drape
- Hand towel

Patient prepping components :

- Forceps
- Sponge stick
- Swabs
- Gallipots

Technical components :

- Dialysis connectors
- Syringes
- Needles
- Blunt fill needles
- NaCl pre-filled syringes or NaCl vials*
- Holder for blood sample
- Clamping forceps
- Tourniquets
- Luer lock caps

Others :

- Dressings
- Adhesive strips
- DASRI bag

*These products are medical devices:

Class IIb, CE marked by MedXL Inc; Notified Body Intertek (CE0413) involved.

Read the labelling instructions carefully.

Class III, CE marked by Becton Dickinson and Company; Notified Body NSAI (CE0050) has been involved.

Read carefully the instructions on the labelling.

Class III, CE marked by CDM Lavoisier; Notified Body GMED (CE0459) has been involved.

Read carefully the instructions on the label.



This composition is an example, we adapt to your protocol.

Dialysis Fistula OFF Kit

Protective and hygiene components :

- Gloves
- Gowns
- Face mask
- Table drape
- Patient drape
- Hand towel

Patient prepping components :

- Forceps
- Swabs
- Gallipots

Technical components :

- Syringes
- Needles
- Blunt fill needles
- NaCl pre-filled syringes or NaCl vials*
- Clamping forceps
- Luer lock caps

Others :

- Dressing
- Adhesive strips
- DASRI bag

*These products are medical devices:

Class IIb, CE marked by MedXL Inc; Notified Body Intertek (CE0413) involved.
Read the labelling instructions carefully.

Class III, CE marked by Becton Dickinson and Company; Notified Body NSAI (CE0050) has been involved.
Read carefully the instructions on the labelling.

Class III, CE marked by CDM Lavoisier; Notified Body GMED (CE0459) has been involved.
Read carefully the instructions on the label.



This composition is an example, we adapt to your protocol.

Dialysis Catheter ON kit

Protective and hygiene components :

- Gloves
- Gowns
- Mob cap
- Face mask
- Table drape
- Patient drape
- Hand towel

Patient prepping components :

- Forceps
- Cotton stick
- Sponge stick
- Swabs
- Gallipots
- Alcohol wipes

Technical components :

- Clamping forceps
- Dialysis connectors
- Syringes
- Needles
- Blunt fill needles
- NaCl pre-filled syringes or NaCl vials*
- Catheter securement devices
- Splitted swabs

Others :

- Dressing
- Scissors
- Adhesive strips
- DASRI bag

*These products are medical devices:

Class IIb, CE marked by MedXL Inc; Notified Body Intertek (CE0413) involved.

Read the labelling instructions carefully.

Class III, CE marked by Becton Dickinson and Company; Notified Body NSAI (CE0050) has been involved.

Read carefully the instructions on the labelling.

Class III, CE marked by CDM Lavoisier; Notified Body GMED (CE0459) has been involved.

Read carefully the instructions on the label.



This composition is an example, we adapt to your protocol.

Dialysis Catheter OFF Kit

Protective and hygiene components :

- Gloves
- Gowns
- Mob cap
- Face mask
- Table drape
- Patient drape
- Hand towel

Patient prepping components :

- Forceps
- Cotton stick
- Swabs
- Gallipots
- Alcohol wipes

Technical components :

- Clamping forceps
- Syringes
- Blunt fill needles
- Needles
- NaCl pre-filled syringes or NaCl vials*
- Citraflow® pre-filled syringes
- Luer lock cap
- Needle free connectors
- Splitted swabs
- Catheter securement device
- Dressing pouch for catheter protection

Others :

- Dressing
- Adhesive strips
- DASRI bag

*These products are medical devices:
Class IIb, CE marked by MedXL Inc; Notified Body Intertek (CE0413) involved.
Read the labelling instructions carefully.
Class III, CE marked by Becton Dickinson and Company; Notified Body NSAI (CE0050) has been involved.
Read carefully the instructions on the labelling.
Class III, CE marked by CDM Lavoisier; Notified Body GMED (CE0459) has been involved.
Read carefully the instructions on the label.



This composition is an example, we adapt to your protocol.

Peritoneal Dialysis kit

Protective and hygiene components :

- Gloves
- Gowns
- Mob cap
- Face mask
- Table drape
- Patient drape
- Hand towel
- Absorbent pad
- Hand gel

Patient prepping components :

- Forceps
- Cotton stick
- Swabs
- Gallipots
- Alcohol wipe

Technical components :

- Clamping forceps
- Dressing
- Adhesive strips
- DASRI bag
- Splitted swabs



Care kit for peritoneal dialysis

Protective and hygiene components :

- Gloves
- Gown
- Mob cap
- Table drape
- Patient drape
- Hand towel
- Hand gel
- Splitted swabs

Patient prepping components :

- Cotton stick
- Swabs
- Alcohol wipe
- Gallipots

Others :

- Dressing
- Adhesive strips



This composition is an example, we adapt to your protocol.



Care kits

- Dressing
- Suturing
- Patient prepping

Dressing replacement kit

Protective and hygiene components :

- Gloves
- Mob cap
- Face mask
- Table drape
- Patient drape
- Hand towel

Patient prepping components :

- Tweezers or Kocher forceps
- Cotton sticks
- Swabs
- Gallipots

Technical components :

- Stylet
- Scissors
- NaCl wound flush

Others :

- Dressing
- Adhesive strips
- DASRI bag



Tracheotomy care kit

Protective and hygiene components :

- Gloves
- Mob cap
- Gown
- Face mask
- Table drape
- Hand towels

Patient prepping components :

- Tweezers
- Cotton sticks
- Swabs
- Gallipots

Technical components :

- Splitted swabs
- Tracheotomy swabs
- NaCl or water vial for cannula cleaning

Others :

- DASRI bag



Central venous access care kit

Protective and hygiene components :

- Gloves
- Mob cap
- Gown
- Face mask
- Table drape
- Patient drape
- Hand towel

Patient prepping components :

- Tweezers or Kocher forceps
- Swabs
- Gallipots

Technical components :

- Splitted swabs
- Scissors
- Catheter securement device
- NaCl pre-filled syringe or NaCl vial*
- Blunt fill needle
- Syringe
- Luer lock cap
- Needle for connector

Others :

- Dressing
- Adhesive strips
- DASRI bag

*These products are medical devices:

Class IIb, CE marked by MedXL Inc; Notified Body Intertek (CE0413) involved.

Read the labelling instructions carefully.

Class III, CE marked by Becton Dickinson and Company; Notified Body NSAI (CE0050) has been involved.

Read carefully the instructions on the labelling.

Class III, CE marked by CDM Lavoisier; Notified Body GMED (CE0459) has been involved.

Read carefully the instructions on the label.



Negative Pressure Therapy (NPT) dressing kit

Protective and hygiene components :

- Gloves
- Table drape
- Hand towel
- Hand gel

Patient prepping components :

- Tweezers or Kocher forceps
- Swabs

Technical components :

- Film dressing
- Scalpel

Others :

- DASRI bag



Patient prepping kit

Protective and hygiene components :

- Gloves
- Face mask
- Hand towel
- Hand gel
- Table drape

Patient prepping components :

- Forceps
- Swabs
- Gallipots
- Cotton sticks
- Sponge sticks
- *Surgical brush*
- Absorbent pad
- Muff
- Skin disinfection pad

Others :

- DASRI bag

Italic text refers to components currently under sourcing (not yet available)



This composition is an example, we adapt to your protocol.

Patient protection kit

Protective components :

- Gloves
- Face mask
- Hand towel
- Table drape
- Gown
- Mob cap
- Shoe cover



This composition is an example, we adapt to your protocol.

Mouth care kit

Protective and hygiene components :

- Gloves
- Mask
- Gown
- Hand towel
- Table drape
- Patient aperture drape

Patient prepping components :

- Forceps
- Swabs
- Gallipots

Technical components :

- Tongue depressor
- Probe
- Spoon excavator
- Cotton stick



This composition is an example, we adapt to your protocol.

Suture fitting kit

Protective and hygiene components :


- Table drape
- Patient aperture drape
- Hand towel

Patient prepping components :

- Forceps
- Swabs
- Gallipots

Technical components :

- Metallic tweezers
- Sharp scissors for dressings
- Blunt scissors
- Needle holder
- Syringe
- Hypodermic needles
- NaCl wound flush
- Stitches (strips)



This composition is an example, we adapt to your protocol.

Suture removal kit

Protective and hygiene components :

- Table drape
- Hand towel

Patient prepping components :

- Forceps
- Swabs

Technical components :

- Metallic tweezers
- *Staple removing forceps*
- Scissors

Italic text refers to components currently under sourcing (not yet available)





Infusion kits

- Centrale venous access
- Peripheral venous access
- Neuraxial access space

CVC fitting kit

Protective and hygiene components :

- Gloves
- Gown
- Mob cap
- Face mask
- Table drape
- Patient drape
- Hand towels

Antisepsis :

- Forceps or sponge stick
- Swabs
- Gallipots
- Hand gel

Technical components :

- Syringe
- Infusion set
- Extension line
- 3 ways stopcpcks
- Needle
- NaCl pre-filled syringe or vial*
- Blunt fill needles
- Luer lock caps
- Needle free connector
- Catheter securement devices

Others :

- Dressing
- Adhesive strips or self adherent wrap
- DASRI bag
- Protective case for connections

*These products are medical devices:

Class IIb, CE marked by MedXL Inc; Notified Body Intertek (CE0413) involved.
Read the labelling instructions carefully.

Class III, CE marked by Becton Dickinson and Company; Notified Body NSAI (CE0050) has been involved.

Read carefully the instructions on the labelling.

Class III, CE marked by CDM Lavoisier; Notified Body GMED (CE0459) has been involved.
Read carefully the instructions on the label.



This composition is an example, we adapt to your protocol.

PVC fitting kit

Protective and hygiene components :

- Gloves
- Gown
- Mob cap
- Face mask
- Table drape
- Patient drape
- Hand towel
- Hand gel

Antisepsis :

- Swabs
- Gallipots
- Skin disinfection pad

Technical components :

- Syringe
- Infusion set
- Extension line
- 3 ways stopcocks
- Needles
- NaCl pre-filled syringe or vial*
- Blunt fill needle
- Luer lock cap
- Needle free connector
- IV Catheter
- Scalp vein cannula
- IV Dressing

Others :

- Dressing
- Tourniquet
- Adhesive strips or self adherent wrap
- DASRI bag

*These products are medical devices:

Class IIb, CE marked by MedXL Inc; Notified Body Intertek (CE0413) involved.

Read the labelling instructions carefully.

Class III, CE marked by Becton Dickinson and Company; Notified Body NSAI (CE0050) has been involved.

Read carefully the instructions on the labelling.

Class III, CE marked by CDM Lavoisier; Notified Body GMED (CE0459) has been involved.

Read carefully the instructions on the label.



This composition is an example, we adapt to your protocol.

Huber Needle withdrawal or CVC disconnection kit

Protective and hygiene components :

- Gloves
- Gown
- Mob cap
- Face mask
- Table drape
- Patient drape
- Hand towel
- Hand gel

Others :

- Dressing
- Adhesive strips
- DASRI bag

Antisepsis :

- Forceps or sponge stick
- Cotton sticks
- Swabs
- Gallipots

Technical components :

- Syringes
- Extension line
- Needles
- NaCl pre-filled syringe or vial*
- Blunt fill needle
- Luer lock cap
- Needle free connector
- Catheter securement device
- Disinfectant cap for needle free connector
- Huber needle withdrawal forceps

*These products are medical devices:
Class IIb, CE marked by MedXL Inc; Notified Body Intertek (CE0413) involved.
Read the labelling instructions carefully.
Class III, CE marked by Becton Dickinson and Company; Notified Body NSAI (CE0050) has been involved.
Read carefully the instructions on the labelling.
Class III, CE marked by CDM Lavoisier; Notified Body GMED (CE0459) has been involved.
Read carefully the instructions on the label.



This composition is an example, we adapt to your protocol.

PVC disconnection kit

Protective and hygiene components :

- Gloves
- Gown
- Mob cap
- Face mask
- Table drape
- Hand towel
- Hand gel

Antisepsis :

- Forceps
- Swabs

Technical components :

- Syringe
- NaCl pre-filled syringe*
- Luer lock cap
- Needle free connector
- Catheter securement devices
- Disinfectant cap for needle free connector

Others:

- Dressing
- Adhesive strips
- DASRI bag

*These products are medical devices:

Class IIb, CE marked by MedXL Inc; Notified Body Intertek (CE0413) involved.
Read the labelling instructions carefully.

Class III, CE marked by Becton Dickinson and Company; Notified Body NSAI (CE0050) has been involved.
Read carefully the instructions on the labelling.

Class III, CE marked by CDM Lavoisier; Notified Body GMED (CE0459) has been involved.
Read carefully the instructions on the label.



This composition is an example, we adapt to your protocol.

Kit for subcutaneous injection

Protective and hygiene components :

- Gloves
- Gown
- Mob cap
- Face mask
- Table drape
- Hand towel
- Hand gel

Antisepsis :

- Swabs
- Skin disinfection pad

Technical components :

- Syringes
- Blunt fill needle
- Needle free adaptor and Transfer Device
- IV catheter
- Scalp vein cannula
- Infusion set
- Extension line
- 3 ways stopcocks
- Luer lock cap

Others :

- Dressing
- Tourniquet
- Adhesive strips
- DASRI bag



This composition is an example, we adapt to your protocol.

Neuraxial Catheter disconnection kit

Protective and hygiene components :

- Gloves
- Gowns
- Mob cap
- Face mask
- Table drape
- Patient drape
- Hand towel

Antisepsis :

- Swabs

Technical components** :

- Syringes
- Blunt fill needle
- NaCl vial or pre-filled syringe*
- Extention line
- 3 ways stopcocks
- Luer lock cap
- Needle free connector
- Catheter securement device

Others :

- Film dressing
- Adhesive strips
- DASRI bag

*These products are medical devices:

Class IIb, CE marked by MedXL Inc; Notified Body Intertek (CE0413) involved.

Read the labelling instructions carefully.

Class III, CE marked by Becton Dickinson and Company; Notified Body NSAI (CE0050) has been involved.

Read carefully the instructions on the labelling.

Class III, CE marked by CDM Lavoisier; Notified Body GMED (CE0459) has been involved.

Read carefully the instructions on the label.

**Following ISO 80369 a new standard connection is used for regional anesthesia: NRRfit, some NRRfit components are already available (syringes and fill needles) please contact customer support for more information



This composition is an example, we adapt to your protocol.

Filling and preparation kit

Protective and hygiene components :

- Gloves
- Gowns
- Mob cap
- Face mask
- Table drape
- Hand towel
- Hand gel

Asepsis :

- Swabs

Technical components :

- Syringes
- NaCl pre-filled syringe or vial*
- Blunt fill needle
- Air intake cannula
- Needle
- Infusion set
- Extension line
- 3 ways stopcocks
- Luer lock cap
- Needle free adaptor and transfer device

Others :

- DASRI bag
- Protective case for connection

*These products are medical devices:

Class IIb, CE marked by MedXL Inc; Notified Body Intertek (CE0413) involved.

Read the labelling instructions carefully.

Class III, CE marked by Becton Dickinson and Company; Notified Body NSAI (CE0050) has been involved.

Read carefully the instructions on the labelling.

Class III, CE marked by CDM Lavoisier; Notified Body GMED (CE0459) has been involved.

Read carefully the instructions on the label.



This composition is an example, we adapt to your protocol.



Central venous access fitting kits

- Implantable chamber
- Catheter
- PICCLine
- MIDLine

Central venous access fitting kit

Protective and hygiene components:

- Gown or surgical gown
- Patient drape with aperture
- Mob cap
- Face mask
- *Surgical brush*
- Hand towel
- Hand gel
- Table drape
- Probe cover
- Gloves

Others :

- Scissors
- Tourniquet
- Ultrasound Transmission gel
- Dressings
- Clamping forceps

Patient prepping components:

- Sponge stick
- Swabs and XRD swabs
- Gallipots
- Forceps
- Sponge forceps

Technical components:

- Blunt fill needle
- Hypodermic needle
- Scalpel
- 3 ways stopcocks
- Needle free connector
- Luer lock cap
- Extension
- Syringes
- NaCL vial* or prefilled syringe
- Catheter securement device
- Needle holder

Italic text refers to components currently under sourcing (not yet available)

*These products are medical devices:

Class IIb, CE marked by MedXL Inc; Notified Body Intertek (CE0413) involved.

Read the labelling instructions carefully.

Class III, CE marked by Becton Dickinson and Company; Notified Body NSAI (CE0050) has been involved.

Read carefully the instructions on the labelling.

Class III, CE marked by CDM Lavoisier; Notified Body GMED (CE0459) has been involved.

Read carefully the instructions on the label.



This composition is an example, we adapt to your protocol.



Catheterisation kits

- Drain
- Probes
- Urinary Catheterisation...

Vesical catheterisation kit

Protective and hygiene components :

- Gown
- Hand towel
- Hand gel
- Gloves
- Face mask
- Table drape
- Splitted or aperture drape

Patient prepping components :

- Sponge sticks
- Kocher forceps
- Tweezers
- Swabs
- Gallipots
- Muff
- Sponge forceps

Technical components :

- Lubricating gel
- Sterile water (pre-filled syringe or vial)
- Syringes
- Blunt fill needle

Others :

- Adhesive strips
- Bag DASRI



This composition is an example, we adapt to your protocol.

Probes and drains fitting kit

Protective and hygiene components :

- Surgical gowns
- Patient drape with aperture
- Mob cap
- Face mask
- Gloves
- *Surgical brush*
- Hand towel
- Table drape

Patient prepping components :

- Sponge stick
- Swabs
- Gallipots
- Forceps
- Sponge forceps

Technical components:

- Needle
- Scalpel
- Syringes (ENFit in case of enteral nutrition)
- Needle holder
- Lubricating gel
- IV catheter
- Sterile water vial

Others :

- Scissors
- Tourniquet
- IV dressing
- Clamping forceps
- DASRI bag

Italic text refers to components currently under sourcing (not yet available)



This composition is an example, we adapt to your protocol.

List of legal manufacturers by component



Reference	Name	Class / N° and name ON	Manufacturer
103 304	10ml NaCl Praxiject® prefilled syringe	I Ib, 0413 Intertek	MedXL Inc.
103 465	20ml NaCl Praxiject® prefilled syringe	I Ib, 0413 Intertek	MedXL Inc.
103 501	10x15 IV dressing	I Ia, 2797 BSI	3M Health Care
103 968	Lubricating gel 2.7g	I Ia, 2797 BSI	Optimum Medical Solution Limited
104 371	Water + 10% glycerin prefilled syringe	I s, 2797 BSI	Optimum Medical Solution Limited
104 544	Neoprene gloves	I Ia, 1370 Bureau Veritas Italia spa	RAYs SPA
102 504	MaxPlus® needle free connector	I Ia, 2797 BSI	Becton Dickinson and Company
100 163	6x8 IV dressing	I s, 2797 BSI	Smith & Nephew Medical Ltd
102 719	10ml NaCl Posiflush® prefilled syringe	III, 0050 NSAI	Becton Dickinson and Company
103 463	Caresite® needle free connector	I Ia, 0123 TÜV Sud	B. Braun Medica Inc.
103 884	Statlock® catheter securement device	I s, 2797 BSI	Becton Dickinson and Company
104 370	Grip-Lock® catheter securement device	I s, 0459 GMED	VYGON
103 432	46.7% Citrafow® prefilled syringe	I Ia, 0413 Intertek	MedXL Inc.
103 109	5ml NaCl plastic vial	I Ia, 0459 GMED	Laboratoire GILBERT
104 594	Swab for tracheotomy	I s, 0123 TÜV Sud	Lohmann & Rauscher International GmGH & Co. KG
104 401	1 way infusion set	I Ia, 1370 Bureau Veritas	DIDACTIC S.A.S
104 490	Safety scalp vein cannula 24G	I Ia, 2797 BSI	Becton Dickinson and Company
103 892	Safety IV catheter 22G	I Ia, 2797 BSI	Becton Dickinson and Company
104 761	Spinal Yale needle 22G- 75mm	III, 0318 agencia espanola de medicamentos y productos sanitarios	Becton Dickinson and Company
104 236	10x10 wrapped NW swabs	I s, 0123 TÜV Sud	Henan Zhende Medical Products Co., Ltd.

List of legal manufacturers by component



Reference	Name	Class / N°and name ON	Manufacturer
104 779	Blue male/Female LL cap	Is, 1370 Bureau Veritas	DIDACTIC S.A.S
104 234	7.5x7.5 wrapped NW swabs	Is, 0123 TÜV Sud	Henan Zhende Medical Products Co., Ltd.
103 543	Red LL cap	Is, 1370 Bureau Veritas	DIDACTIC S.A.S
103 992	Needle free connector with Clave®	Ila, 0050 NSAI	ICU Medical, Inc.
104 822	3 ways infusion set	Ila, 0123 TÜV Sud	Jiangsu Webest Medical Product Co., Ltd.
104 820	Y infusion set with flow controller	Is, 0123 TÜV Sud	Sendal S.L.
104 778	I3 Y infusion set with Clave®	Ila, 0050 NSAI	ICU Medical, Inc.
104 840	50 ml NaCl plastoc vial	Ila, 0459 GMED	Laboratoire GILBERT
104 841	Sharp container 3l- pas dans les sets	Non DM	Frontier Plastics Limited
104 859	ANIOSAFE® soap pod	Non DM	ANIOS
104 845	CHEMOCLAVE® Universal vented vial spike	Is, 0050 NSAI	ICU Medical, Inc.
104 867	Isolation gown size L	I	DISPOMED
104 868	Isolation gown size XL	I	DISPOMED
104 870	Self-adherent wrap (Coban®)	I	3M Company
104 864	Eclipse® needle 21G	Ila, 0050 NSAI	Becton Dickinson and Company
104 865	Eclipse® needle 25G	Ila, 0050 NSAI	Becton Dickinson and Company
104 918	3 ways stopcock with 30 cm extension line	Ila, 0123 TUV SUD Produt Service GmbH	JIANGSHU WEBEST
104 417	Safety IV catheter 22G with Blood control	Ila, 2797 BSI	Becton Dickinson Infusion

This image shows a full page of white paper with horizontal dotted lines. The lines are evenly spaced and run across the width of the page, providing a guide for handwriting practice. There are no margins, text, or other markings on the page.

About Hemodia

Founded in 1985, Hemodia has developed expertise in medical devices. Based on core business, medical tubing and polymer components, Hemodia has rapidly become a major player on the custom-made care kits market. From dialysis to anesthesia sets, Hemodia adjusts its offer to every request and speciality : oncology, infusion, hygiene...

Thanks to the concept of internalisation of conception, production, sterilisation, storage and distribution – Hemodia anticipates to the customers' needs and answers to all requirements concerning deliveries once orders have been placed. With production in France and Tunisia it ensures total control on supplies and equipment.

Thanks to quality management systems, Hemodia stands by your side by reducing the risk of infections, improving compliance and best practice, avoiding wastage and creating value. The high level of requirements guarantees both professionals and patients the best quality products.

Trusting in Hemodia is choosing knowledge and recognised expertise.



For more information, please contact :

Customer service

+33 (0)5 61 00 71 81
hemodia@hemodia.com

www.hemodia.com



Hemodia SAS

104 rue du Chêne vert
31670 Labège - France

